

Etac Stapel bathroom stool

Stapel – a reliable bathroom stool

Stapel is a sturdy and safe stool which can be used in many situations. It is ideal to use when stepping in and out of the bathtub or as a footstool when using the toilet. Stapel has a plastic top that easily can be snapped off for cleaning.



Stapel makes stepping in and out of the bathtub safe.

Accessory



Telescopic leg

83010002
5–10.5 cm (2–4¼"), 4 pcs
Item no. 83010001
5–18 cm (2–7"), 4 pcs

Description

Stapel 81204002

Dimensions

Height: 18 cm (7")
Footprint: 48.5 – 49 cm (19 – 19¼")
Standing area: 38 – 43 cm (15 – 17")

Material

Frame: Powder-coated tubular steel
Seat: Polypropylene

Weight

3 kg (6.6 lbs)



Clean the product with a solvent-free cleaning agent (pH level 6–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.