

Etac Swift Mobile 24"


A perfect solution for independency.

Facilitate side transfer

Swift Mobile 24 is easy to use and propel. The rear wheels are perfectly placed to be easy to reach and give the chair a high manoeuvrability. Thanks to the seat design there is sufficient space in front of the wheels to make a safe side transfer. Furthermore, the seat height can be set in an appropriate height that facilitate the transfers.



A standard Swift Mobile can also be remodeled to a Swift Mobile 24" by assembling the accessory "Rear Wheel Kit". See page 51.



Accessories

See our wide range of accessories on page 50-51.

Description

Swift Mobile 24" assembled 80209402
 Swift Mobile 24" assembled, incl. pan holder 80209408
 Swift Mobile 24" incl. bucket holder 80209407
 Standard: back support, seat, arm supports, foot supports and 4 lockable castors, rear wheels 24",brakes.

Material

Chassis: powder painted steel
 Seat/Arm supports: polypropylene
 Back support: polyester/Velcro
 Foot support: polypropene
 Castors: polyamid/TPE, lockable, ball bearing
 Rear wheels: solid tyres, polyamid, ball bearing
 Brakes: polyamid, stainless steel

Dimensions

Total width: 74 cm (29")
 Total length: 106 cm (42")
 Adjustable seat height: 45–60 cm (17 3/4–23 1/2")
 Width between arm supports: 48 cm (19")

Weight

Weight: 24.4 kg (53 lbs)



Clean the product with a solvent-free cleaning agent (pH level 6–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.

Accessories for all Swift Mobile models



Head Support
80209447
Complete with attachment. Adjustable in height, depth and sideways. Material: Polyuretane, stainless steel and zinc. max. 85° C (185°F) or with alcohol wipe



Side support
Grey: 80209509
Supports the upper trunk.
Material: Polyester fabric, Velcro, polyethylene.
Washing at max 60°C (140°F)



Arm supports widening kit
Pair: 80209441
Increases the space between the arm supports by 6 cm (2½"). Standard for Swift Mobile/Tilt 160.



Cross bar
Length: 60.5 cm (24.2") 80303034
Length: 66.5 cm (26.6") 80209452
For safety and additional support.
Material: Aluminium, polypropene and cellular polyethylene. Notice: 66.5 cm is for Swift Mobile 160 and Swift Mobile Tilt 160.



Hemi Arm Support
Pair: 80209035
Fits both left and right side. Material: Cellular polyethylene, aluminium, stainless steel, polystyrene.



Chest harness strap
80209442
Adjustable in length. Attach to back support. Always use together with positioning belt.
Material: Polyester



Arm supports lock
Pair: 80209445
Locks the arm supports in down-folded position. Material: Stainless steel.



Positioning belt
80209440
Adjustable in length. Attachment with velcro straps.
Material: Polyester



Comfort seat
80209433
A seat pad which reduces the recess opening to 18 cm (7"). Material: Foamed polyurethane. max. 85° C (185°F) or with alcohol wipe



Tilter
80209449
For easy transfer over thresholds and other obstacles



Seat cover
80209434
When a complete seat is required. Material: Foamed polyurethane. max. 85° C (185°F) or with alcohol wipe



Heel straps
Pair: 80209432
Adjustable in length.
Material: Polyester



Soft Back Pad
80209448
Warm, soft and water repellent. Material: Polyethylene. max. 60° C (140°F) or with alcohol wipe



Elevating Leg- and calf support
80209437
Fits both left and right side.
Material: ABS



Splash Guard
80209430
Fits into the recess opening. To be placed when user sits in the chair.
Material: Polyurethane




Foot support for leg- and calf support
Pair: 80209443
Fits both left and right side. Foot support adjustable in 6 positions.
Material: Stainless steel, polypropylene.



The comfortable foot supports provide stability and protects the feet from the castors.




Attachment pan- and bucket holder,
Grey, set: 80209435
For pan- and bucket holder. To be mounted on the lower side of the seat.
Material: Polyamide.



Soft cushion with standard opening
80209065
Material: Polyether/Sealed coated cover



Pan holder
Grey, pair: 80209257
Fits both Swift Mobile pan and other with thin rim. Order attachment to pan holder separately. Material: Polypropene.



Soft cushion with oval opening
80209067
Material: Polyether/Sealed coated cover



Pan
Grey: 80209254

Pan, incl lid with handle
Grey: 80209255

Lid with handle
Grey: 80209256
Material: Polyamide.



Anti tip, pair
Excluding attachment (for Swift mobile 24"): 80209438
Including attachment (for Swift Mobile): 80209439
(for Swift Mobile Tilt): 80209459
Increases safety if the center of gravity is far back, e.g. leg amputations.
Material: Aluminium



Attachment for low mounting of pan/bucket
Grey, set of four: 80209259
Material: polypropylene



Rear wheel kit
80209436
24" rear wheels with solid tires. With built in attachment for anti tip.
Material: Polyamide, stainless steel.
Incl. attachment and brakes. Do not fit Swift Mobile 160 or Swift Mobile Tilt.



Leg support
80209446
Leg support with foldable and angle adjustable footplate.



Power tilt kit
80209499
For Swift Mobile Tilt and Swift Mobile Tilt 160 only. Hand control and charger are included.