

Etac Swift Commode

A three in one product with many smart functions.

Comfortable and functional

Swift Commode is comfortable and discreet in the bedroom. Under the seat padding hides a removable pan with a tightly securable lid that locks automatically when lifted. Swift Commode can be set with a slight forward tilt to facilitate easier access. The soft seat and back support paddings are standard.

3 in 1

Swift commode can also be used as a shower chair or as a freestanding toilet seat raiser. The soft ferrules give an excellent grip on all surfaces.

Easy cleaning

Swift Commode is designed with easy cleaning in mind. The smooth surfaces, the easy detachable paddings and the fact it is corrosion free are important features to make cleaning easier. Swift Commode comes flat packed and is easy to assemble and height adjust, without using tools.



With the paddings, the Swift commode has a discrete design and is comfortable also for dressing and undressing.

Description

Swift commode 81702030

Dimensions

Total width: 56 cm (22")
Seat width: 54 cm (21¼")
Seat height min–max: 42–57 cm (16½–22½")
Width between the arm supports: 45 cm (17¾")

Material

Seat, back and arm supports: polypropylene
Legs: aluminum and polyamide
Ferrules: TPE
Back and seat pads: Polythylene with nylon fabric

Weight

5.4 kg (12 lbs)



Clean the product with a solvent-free cleaning agent (pH level 6–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.



Integrated pan lid

The pan with lid fits smoothly and becomes one unit with the seat.

Accessory



Soap dish

80209266
Handy for toilet paper
Material: polypropylene