

Etac Molift Mover 180

Etac Molift Mover 180 is a compact and Light weight lifter, suitable both in home care and institutional environments.



An everyday favourite

The Molift Mover 180 has a lightweight alloy construction that makes the lifter very easy to maneuver. The large lifting interval, 27 – 163 cm (10.6" – 64.2") allows lifting from high surfaces as well as lifting from the floor. Molift Mover 180 has also the option of 2 – or 4-point suspension.

Ergonomical handling

For improved carer ergonomics, the Molift Mover 180 is equipped with Molift's unique push bar that provides multiple and flexible grip possibilities. It is also equipped with Molift's environmentally friendly NiMH-battery.

Molift Mover 180, like all Molift lifters, is also equipped with the unique service notification system that notifies when service is required. The integrated service software can also be used to retrieve data about lifter usage and lifting patterns.

Features

Integrated service software

Automatic service notification system determines when the lifter requires service and maintenance.



Unique push bar

The unique push bar provides the best natural grip, creating good ergonomics and many grip possibilities.



Many options

Mover 180 can be used with both a 2-point and a 4-point suspension. They come in six variants and can be replaced easily without tools.



Easy to maneuver

Molift Mover 180 has a total weight of only 32 kg (70.5 lbs) which, together with the design, makes it easy to maneuver.



You are here: [Transfer solutions](#) / [Lifts and Mobile Lifts](#) / [Molift Mover 180](#)
