

Etac Molift Mover 205

All-round mobile lift for all situations.



Find optional suspensions and detailed information at page 26–27.

A workhorse for frequent and everyday use

Etac Molift Mover 205 is a flexible all-round mobile passive lift for all lifting and transfer situations in nursing, institutional and hospital settings. Although having a low total weight of 39.3 kg (89.5 lbs) facilitating wheeling and maneuverability, the unique construction allows a safe working load 205 kg (452 lbs).

Excellent reach, low and high

Molift Mover 205 has an excellent lifting height, with the lifting interval 37 – 179 cm (14.5" – 70.4"). It has a standard 4-point suspension that used together with a Molift Sling provides a comfortable and spacious position for the user. A quick release function provides options of other suspensions in different dimensions. Equipped with environmentally friendly NiMH-battery.

Description

Etac Molift Mover 205 excl. suspension
Standard 14000
Low base 14050
Support arms 14150

Material

Aluminium, steel

Service Software

Included

Total weight

Regular: 39.3 kg (89.5 lbs)
Low base: 42,6 kg (94 lbs)
With support arms: 45.6 kg (100.5 lbs)

Height

139 cm (54.9")

Length

122 cm (48.0")

Width chassis/legs outer

70 cm (27.6")

Width open legs inner

107 cm (42.1")

Lifting height

37 – 179 cm
(14.5" – 70.4")

Lifting interval

146 cm (57.5")

Leg height

11.5 cm (4.5")
Low base: 6 cm (2.4")
With support arms: 11.5 cm (4.6")

Battery

NiMH 26.4 V – 2.2 Ah

Battery charge time

Ca 3 hours

Lifting speed

4 cm/ sec
(1.6"/ sec)

Rated performance

50 lifts at: 75 kg and 50cm
(165 lbs and 19.7")

Cleaning instructions

Use pH-neutral cleaning agents. If disinfection is required, use isopropyl alcohol. Do not use abrasive or corrosive substances for cleaning.



Features



Walking arms

Molift Mover 205 is also available with additional Support Arms. This in combination with the excellent maximum height provides an ideal solution for gait training.



Extra low base

Etac Molift Mover 205 is also available with an extra low base chassis. A base height of only 6 cm is one of the lowest in the market.



Many options

Mover 205 can be used with both a 2-point and a 4-point suspension. They come in six variants and can be replaced easily without tools.



Unique push bar

The unique push bar provides the best possible natural grip, creating good ergonomics and many grip possibilities.



The Etac Molift Stretcher can be used directly on the standard 4-point suspension, making it easy to lift a multiple trauma patient, a patient with suspected spinal cord injuries, or an unconscious patient.