

There is Power Mobility. And then there's

POWERFUL MOBILITY



Spring



There is Power Mobility. And then there's

POWERFUL MOBILITY

Spring

The Heartway Spring is a manual wheelchair with a lightweight aluminum frame. The Spring uses a Tilt N Space seat that operates manually. The Tilt N Space seat reclines up to 30 degrees.

The Spring Has a weight capacity of 300 lbs and weighs only 84 lbs. The Spring features 24" quick-release drive wheels. The ergonomic seat is available in 16", 18", and 22" widths.



MAX SPEED
N/A



LOAD CAPACITY
300
LBS



BATTERY
N/A



RANGE
N/A

Length	45"
Width	26" (16" seat width) 27" (18" seat width) 31" (22" seat width)
Seatback Height	22" (16" & 18" seat widths) 24" (22" seat width)
Overall Height - at Minimum Tiller Height	N/A
Legroom - Floorboard Length x Width	16-20"
Deck Height - Floorboard to Floor	4-8"
Turning Radius	39"
Weight	84 lbs

Seat Type	Rehab Seat
Seat - Width & Depth & Height	16" / 18" / 22"
Seat - Cushion Thickness	4"
Seat Height - Cushion to Deck	16"
Seat Height - Cushion to Floor	22"
Tire Type	Solid Tires
Rim Type Rear	Industrial Plastic
Rim Type Front	Industrial Plastic
Outside Diameter Rear	24" in

Tire Width Rear	1-3/8"
Outside Diameter Front	8"
Tire Width Front	2"
Climbing Angle - Maximum Recommended	6 degree
Load Capacity	300 lbs
Colors	Silver
Elevating Leg Rests	Standard

(*) SPECIFICATION PERFORMANCE IS DETERMINED IN THE LABORATORY SETTING BASED ON IDEAL CONDITIONS (MINIMAL FRICTION, NO WIND AND CONTROLLED TEMPERATURE AND PRESSURE WITH AN AVERAGE STANDARD WEIGHT). ACTUAL PERFORMANCE WILL BE DIFFERENT THAN RESULTS OBTAINED DURING LABORATORY TESTING