



TRANSITIONS® Angled Entry Mat

Technical Specifications

HEIGHT & USAGE GUIDELINES (RISE)	1.5"	2.5"
MADE FROM 100% RECYCLED RUBBER	✓	✓
SLIP-RESISTANT SURFACE	✓	✓
AVAILABLE COLORS	Black or Grey	Black or Grey
USABLE SIZE (H x L x W)	1.5" x 14" x 36"	2.5" x 25" x 36"
OVERALL SIZE (H x L x W)	1.5" x 14" x 40"	2.5" x 25" x 43"
WEIGHT	14 lbs.	38 lbs.
WEIGHT CAPACITY (DO NOT EXCEED)	850 lbs.	850 lbs.
WARRANTY	Lifetime Limited	Lifetime Limited
COUNTRY OF ORIGIN	China	China

- Weights and dimensions are approximate and subject to change.

- If using with mobility equipment, refer to your mobility equipment's owner guide for acceptable usage, including proper slope/incline, chair direction, etc.; never exceed its recommendations.

