



# Freedom to move.

**Opus**<sup>360</sup>

**NASAL PILLOWS MASK**

The Opus<sup>360</sup> Nasal Pillows Mask is small, quiet and superbly comfortable, offering unprecedented freedom. The slim-line design allows for almost any sleeping position without mask displacement. Its innovative 360-degree pivoting ball-and-socket elbow provides greater versatility. With its new Silicone Pillows and facially Contoured Mask Frame for extra stability, patient comfort is paramount.


Elegant in design and simplicity, let the Opus<sup>360</sup> Nasal Pillows Mask from Fisher & Paykel Healthcare work for you today!

**Contoured Mask Frame  
& "Off-the-Lip" Pillows**



**Fisher & Paykel**  
HEALTHCARE

# Direct Nasal Delivery – Opus<sup>360</sup>



The Opus<sup>360</sup> Nasal Pillows Mask is small, quiet and superbly comfortable, with anatomical pillow-style Silicone Seals and a facially Contoured Mask Frame. A ball-and-socket elbow combines with an innovative anchoring system to give unprecedented freedom of movement and stability. The Opus<sup>360</sup> Nasal Pillows Mask is a complete, Out-of-Box Solution, fully assembled and ready for use with alternative sizes of Silicone Pillows included.

## UNIQUE MASK FEATURES FOR GREATER FREEDOM, COMFORT AND COMPLIANCE



### Anatomically Shaped Pillows Mask comfort

The Anatomically Shaped Pillows provide easy fitting with an excellent seal. The contours of the nasal pillows ensure greater comfort. The range of sizes provided (S, M, L) ensures comfort and seal are not compromised.




### Contoured Mask Frame Mask comfort and stability

The facially Contoured Mask Frame, with smooth soft pads, conforms naturally to the shape of the face to provide a comfortable fit and reliable support. Its slim-line design allows side-sleeping without mask displacement.




### Tube Anchoring System Mask stability

The innovative pivot system allows the tube to be anchored at the top or side, or left free to track. Elegant in design and simplicity, it firmly holds the tube in the desired position. Alternatively, the tube can be left unanchored so that the ball and socket elbow allows true 360-degree freedom.



### Efficient Exhaust Quiet usage

The advanced air diffuser system is quiet, minimizing the noise and draft of the exhaust flow to prevent both sleeper and partner disturbance.



### "Off-the-lip" Pillows Mask comfort

The Silicone Seal Pillows are soft on the nose and designed to be "off the upper lip" for maximum comfort.



### Premium Frosted Silicone Seal Optimizing seal technology

The enhanced seal technology has a soft, smooth feel and slides easily into position. Its frosted texture optimizes mask seal performance.



### Hidden Strap Release One-step removal

The easy one-step removal of headgear for cleaning enables reassembly without having to make readjustments. Simply slide the Side Strap down and out of the Hidden Strap Release hooks.

## CONTENTS

- Complete Mask with Mask Frame assembled and ready for use
- All three sizes of Silicone Seal included (S, M, L)
- User-friendly Fitting Guide

REORDER CODE: HC482