

# Find Your Freedom!

ACCEPTANCE • INDEPENDENCE • DISCOVERY • BRAVERY  
ACCOMPLISHMENT • BALANCE • DRIVE • FRIENDSHIP • AGILITY • BRAVERY  
PERSEVERANCE • DETERMINATION • RESILIENT • FULLFILLMENT •  
COURAGE • BELIEVING • VALOR • WILLPOWER • GOAL-SETTING  
BRAVERY • AGILITY • RECOGNITION • DETERMINATION • BALANCE • RESILIENCY • B  
WILLPOWER • GOAL-SETTING • BALANCE • TRIUMPH • SELF RELIANT • MOTOR SKILLS

# **FREEDOM** *Concepts Inc.*



# DISCOVERY SERIES



## Discovery is the Mobility Device that started it all!

Today, it's the bike that children and young adults want for its funky design, bright colors, and great ride. More than just freedom, Discovery lets individuals create their own sense of cool.

From those with simple balance problems, to those more involved cases of Cerebral Palsy, Spina Bifida, Down syndrome, Head Injury, Muscular Dystrophy, Autism, and Dwarfism—the Discovery Series offers the perfect combination of therapy and fun.

With quick release mechanisms, the Discovery adapts easily from one rider to another, making it the first choice of many therapists, schools, and care facilities. The Discovery is designed to provide excellent stability and to fit through standard door openings.

The lightweight, agile Discovery Series rides on 12" and 16" wheels.

### DISCOVERY SERIES

Rider Age .....	18 mon.- 14 years
Rider Inseam....	6" - 26"
Max Load .....	60 lbs - 150 lbs
Bike Weight.....	40 lbs - 45 lbs
Wheel Size .....	12" - 16"
Length .....	42" - 48"
Width .....	29"



## Available With Lots Of Standard Features:

- Adjustable handlebars can move completely out of the way for dismounting.
- Adjustable seat in low back or high back styles comes complete with safety belts.
- A unique bottom bracket and crank can be adjusted in all directions to accommodate different hip to feet positions.
- A discreet **Rear-Steer Assistance Handle** lets the rider learn to ride and stop with a little help.
- Direct drive and low gear ratios make pedaling easy for a beginner, yet challenging for an aggressive rider.
- Velcro straps keep feet on the pedals or supporting Neoprene or ABS footplates.
- Locking hand brake mechanisms make for easy mounting and dismounting.
- Powder-coated color options.

## Additional Options To Consider:

- For additional stopping power, rear-brake controls can be mounted on the handlebars, or on the **Rear Steer Assistance Handle™** for caregiver operation.
- Velcro gloves provide the rider superior hand control.
- **TEC Drive™**—offering easy, average or challenging pedalling power.
- Spoke covers, enclosed chain guards and parking brakes are just some of the safety features to consider when ordering your F3 Bike. *Be sure to use your orange safety flag included with your purchase.*

### DCP MINI

Frame color shown: Rainbow with NEW Jungle Decal Package. Available with: An enclosed chain guard. Constructed of ABS plastic, the chain is enclosed to protect extremities from the chain and sprocket.



### DCP 12

Frame color shown: Hot Rod Purple with NEW Beach Decal Package. Available with: Boot style Neoprene Adjustable or Molded ABS footplates to aid in proper positioning.

(See Product Size/Features Guide for details.)



### DCP 16

Frame color shown: Candy Blue. Available with: Solid (*never go flat*) tires. Never pump a tire with air again. Great for multi-user environments as a maintenance free, no-scoff alternative.

(Tire sizes 12" and 16")

**“The thrill of seeing a physically disabled child sit safely on a bicycle is priceless.**

To watch legs that really don't have the ability to pedal actually rotate around on the Freedom Concepts bike is extraordinary. Thanks for making this little six year old's face light up as he feels like one of the crowd.” — William's Mom, Heidi



# Discovery Series: Rider Spotlight



## No Bike Like It

Freedom Concepts bikes have affected us in a huge way. It's kept Waverly in really good shape, it's kept her strength, and I really feel it's keeping her strong.

*I don't think there is any other brand of bike that would work.*

She had tried some bikes that were kind of "supped-up," they were the regular bikes that were made to be a little more stable, but she wasn't able to get it peddling herself.

Waverly has a seizure disorder, her bike is made so it has straps around the chest, so incase she does have a seizure or she's having a weak day, she's not going to fall off the bike, she's not going to injure herself because the bikes so stable,

*nowhere else we'd be able to get a bike like that.*

- Colleen Leduc (Waverly's Mother)

Discovery Mini is available with optional **Rear Platform.** Made of durable ABS Plastic.



For those with allergy concerns, all F3 products are **Latex free.**



### DCP 16 WG/EX

4" wider than our standard Discovery bike for added stability. Also available in an EX frame that is extended 9" longer. Comes with 12" or 16" tires.

*(Recommended for individuals with Dwarfism.)*



### TEC Bike™

Features the TEC Drive™ 3-speed reciprocator allowing for on-the-fly gear change ability from easy, to average, to challenging pedal power.

*(Shown with the JoggUR™ Handle Rear-Steer Assistance.)*

