

P  
O  
W  
E  
R  
C  
H  
A  
I  
R



Seat Colors  
Charcoal Gray

Body Colors  
Candy Apple Red

# Roll

**New!**

## <sup>GOLDEN</sup> *Alanté* **SPORT**

The Alanté name continues to mean quality and value! The first thing that will catch your eye with the all-new Alanté Sport is the unique black and red shroud design that is stylish and functional. Featuring a full size, captain's seat, the Alanté Sport is a front wheel drive power base with our black, non-marking drive tires on high strength steel rims and shiny hubcaps.



Rise. Roll. Rest. Ride.



**YOUR LIFE IN MOTION**

## Alanté Sport GP-205 Specifications

<b>Model Number</b>	<b>GP205CC</b>
Medicare Code	K0823
Weight Capacity	300 lbs
Drive Wheels	Front
Maximum Speed	3.5 mph
Operating Range <sup>1</sup>	16.5 miles
Ground Clearance	1.75"
Turning Radius	25"
Type Batteries	2-U1
Colors	Red
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	PG VR-2 50 Amp Remote
<b>Dimensions:</b>	
Length	40.5"
Width	22.25"
Height (ground to top of back of seat)	37" - 39"
Height (ground to top of back of headrest)	44.5" - 48"
Ground to Top of Seat	20.5" - 22.5"
Number of Seat Height Adjustments	3
Size of Increments	1"
<b>Weight of Unit: (assembled)</b>	
Base	67 lbs
Seat w/Arms	46 lbs
Batteries	22 lbs ea.
<b>Tires:</b>	
Front	9"
Casters	6"
Front Anti-Tip Wheels	4"
<b>Standard Seat:</b>	
Back Height - no headrest	19"
Back Height with headrest	26.5" - 30"
Width x Depth	18" X 18"
Color	Gray Vinyl
<b>Optional Seat:</b>	
Back Height - no headrest	21"
Back Height with headrest	28.5" - 32.5"
Width x Depth	20" x 18"
Color	Coffee/Sand, Charcoal/Gray
<b>Footrest</b>	
Height adjustable	Yes
Angle adjustable	Yes
<b>Warranty:</b>	
Frame	5 Years
Drive Train	13 Months
Electronics <sup>2</sup>	13 Months

<sup>1</sup>. Battery range will vary due to rider weight, drive surface and terrain.

<sup>2</sup>. Electronics warranty excludes batteries.

## Features

- 25" turning radius
- 300 lb weight capacity
- 9" driving wheels with alloy rims
- Enlarged 4" anti-tipper wheels

**5 Year Warranty:**  
Power Chair Frame

**13 Month Limited Warranty:**  
Electronic Controller  
Drive Train Components

## Power Chair Accessories and Options



Pack-N-Go



Quad Cane Holder




Crutch Holder



Rear Basket



Oxygen Holder



Walker Holder



Cup Holder



Disassembles easily into 6 pieces!