

LIGHTWEIGHT PORTABLE SCOOTER



Luxury Meets Rugged in TWO NEW Powerful Scooters!

USB Charging Port Standard!

Seat Color
Black Vinyl

Body Colors
Blue
Red



3 wheel GB119
4 wheel GB149

Introducing America's first full-time luxury travel scooter! Enjoy upgraded features in a whole new class of scooters in our Buzzaround Series! Full-size comfort comes with the convenience of true portability on the Buzzaround LX, along with standard features that include a stylish captain's seat, an infinite adjustable tiller at your fingertips, and a USB charging port right under the dash board. With a 375 lb. weight capacity, comfort suspension springs, 40.5" turning radius, and up to 18 mile range on two U1 batteries, the Buzzaround LX offers ample foot room to take you where you want to go in comfort and style!





Full-size comfort with the convenience of true portability!

| Item | GB119 | GB149 |
|------------------------------|-----------|-----------|
| Medicare Code | Not Coded | Not Coded |
| Weight Capacity | 375 lbs. | 375 lbs. |
| Drive Wheels | Rear | Rear |
| Maximum Speed | 5 mph | 5 mph |
| Operating Range ¹ | 18 miles | 18 miles |
| Ground Clearance Under Deck | 4" | 4" |
| Turning Radius | 40.5" | 54" |
| Type Batteries | 2-U1 | 2-U1 |
| Colors | Red, Blue | Red, Blue |
| Freewheel Mode | Yes | Yes |
| Electronic Speed Control | Yes | Yes |
| Electro-Mechanical Brakes | Yes | Yes |
| Charger | Off-board | Off-board |
| Controller | 90 Amp | 90 Amp |

Dimensions

| | | |
|-------------------------------------|---------------|---------------|
| Length | 45" | 44" |
| Width | 22" | 22" |
| Height (ground to top of seat back) | 37"-39" | 37"-39" |
| Height (ground to top of headrest) | 43.5"-48.5" | 43.5"-48.5" |
| Ground to Top of Seat | 21.75"-23.75" | 21.75"-23.75" |
| Ground to Top of Deck | 6.75" | 6.75" |
| Deck to Top of Seat | 15.75"-17.75" | 15.75"-17.75" |
| Number of Seat Height Positions | 3 | 3 |
| Size of Increments | 1" | 1" |
| Front Axle to Rear Axle | 34.7" | 33.8" |

Weight of Unit:

| (assembled) | 168 lbs. | 174 lbs. |
|-------------|-------------------------------|-------------------------------|
| Front Half | 50 lbs. | 57 lbs. |
| Rear Half | 36 lbs. | 36 lbs. |
| Seat w/Arms | 32 lbs. | 32 lbs. |
| Batteries | 24 lbs./pack x 2 = 48 lbs. | 24 lbs./pack x 2 = 48 lbs. |

Tires:

| | | |
|----------------------|--------------------|--------------------|
| Standard Tires | Solid, low profile | Solid, low profile |
| Optional Tires | N/A | N/A |
| Front | 9" | 9" |
| Rear | 9" | 9" |
| Rear Anti-Tip Wheels | 2" | 2" |

Standard Seat

| | | |
|----------------------------------|-------------|-------------|
| Seat Back Height - no headrest | 17" | 17" |
| Seat Back Height - with headrest | 23"-26" | 23"-26" |
| Width x Depth | 17" x 16" | 17" x 16" |
| Color | Black Vinyl | Black Vinyl |

Optional Seat:

| | | |
|---------------------------|-------------|-------------|
| Back Height - no headrest | 16" | 16" |
| Width x Depth | 20" x 18" | 20" x 18" |
| Color | Black Vinyl | Black Vinyl |

Warranty:

| | | |
|--------------------------|-------------------|-------------------|
| Frame | Lifetime Warranty | Lifetime Warranty |
| Drive Train | 2 Year Warranty | 2 Year Warranty |
| Electronics ² | 2 Year Warranty | 2 Year Warranty |
| Batteries | 1 Year Warranty | 1 Year Warranty |

¹ Battery range at 200 lbs., and will vary due to rider weight, drive surface, terrain, and battery type.

² Electronics warranty excludes batteries.



1. An Infinite Adjustable Tiller (no knuckle) at your fingertips!



2. The tiller folds easily and locks in place.



3. A stylish new profile offers ample leg room!



4. Ultrabrite LED headlight and taillight!



5. USB & Battery Charging Ports conveniently located on the tiller under the control panel.

Features

- Infinite Adjustable Tiller
- USB Charging Port
- Stylish Captain's Seat
- 375 lb. Weight Capacity
- Comfort Suspension Springs
- 40.5" Turning Radius
- Up to 18 Mile Battery Range



Disassembles into 6 pieces!
3 wheel model shown