



Seat Colors
Black Vinyl

Body Color
Blue
Red

**Disassembles
Quickly &
Easily!**

LiteRider^{GOLDEN}
ENVY

The all new LiteRider Envy Power Wheelchair is both stylish and durable with a tubular design that is unique to the Golden LiteRider family! The LiteRider Envy, featuring 2-22 amp hour batteries for an operating range of up to 15.5 miles, is Golden's Retail Cash Power Wheelchair! The LiteRider Envy has an impressive turning radius of 26", and like the other LiteRiders in the family, it disassembles quickly and easily for transport. The heaviest piece is only 35 lbs., making storage convenient in the trunk of any vehicle. When you are looking for a compact power wheel chair that is fun to drive and easy to transport – look no further than the Golden LiteRider Envy!

 **GOLDEN**

Model Number	GP162
Medicare Code	Not Coded
Weight Capacity	300 lbs.
Drive Wheels	Rear
Maximum Speed	3.5 mph
Operating Range ¹	15.5 miles
Ground Clearance	2.5"
Turning Radius	26"
Type Batteries	2-22 AH
Colors	Red, Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	PG VR-2 50 Amp Remote

Dimensions:

Length	33"
Width	22.5"
Height (ground to top of back of seat)	33" - 34"
Height (ground to top of back of headrest)	N/A
Ground to Top of Seat	20.5" - 21.5"
Deck to Top of Seat	18" - 20"
Number of Seat Height Positions	2
Size of Increments	1.0"

Weight of Unit: (assembled)

120 lbs.	
Front Half	34 lbs.
Rear Half	35 lbs.
Seat w/Arms	22 lbs.
Batteries	29 lbs/pack

Tires:

Flat Free	
Drive Wheels	9"
Casters	6"
Anti-Tip Wheels	3"

Standard Seat:

Stadium Style	
Slide Seat	N/A
Back Height - no headrest	16"
Back Height with headrest	N/A
Width x Depth	17" x 16"
Color	Black Vinyl

Footrest:

Yes	
Height adjustable	Yes
Angle adjustable	Yes

Warranty:

Five Years	
Frame	Five Years
Drive Train	13 Month Warranty
Electronics ²	13 Month Warranty

¹ Battery range tested at 200 lbs. and will vary due to rider weight, drive surface and terrain.

² Electronics warranty excludes batteries.

Disassembles easily for transport!



1. Remove seat



2. Pull out battery pack



3. Grasp the frame and locking handle; lift the handle to unlock, and lift the frame up and off the drive train.



Perfect for transport!

Available Accessories:

- Cane Holder Tube
- Safety Flag
- Oxygen Tank Holder
- Power Chair Cover
- Cup Holder
- Touch-up Paint
- Quad Cane Holder
- Deluxe Pack N' Go
- Walker Holder
- Seat Belt (Optional)
- Crutch Holder

Some accessories require a mounting bracket or mounting clips.