



CHERRY RED



PARADISE BLUE



SOLAR FLARE ORANGE  
*Limited Edition*



ENVY GREEN  
*Limited Edition*



SATIN SILVER  
*Limited Edition*



SUNBURST YELLOW  
*Limited Edition*

**Seat Color**



Black Vinyl



**Dynamic LiNX Controller  
Now Standard, Plus a New  
Premium Seat Option!**

*Available in 6 Vibrant Colors  
with a Larger Foot Plate  
and New  
Shroud Design*

**Disassembles  
Quickly &  
Easily!**

**LiteRider<sup>GOLDEN</sup>  
ENVY**

**GP162**


The all-new LiteRider Envy Power Wheelchair is both stylish and durable with a tubular design that is unique to the Golden LiteRider family! A Dynamic LiNX controller comes standard and the traditional style makes driving simple and reliable. The new premium seat offers contoured support with an increased back height of 17" for a longer seating tolerance and added comfort. The Envy, featuring 2-22 amp hour batteries for an operating range of up to 15.5 miles, is Golden's Retail Cash Power Wheelchair! It has a larger footplate, an impressive turning radius of 28.5", and disassembles quickly and easily for transport. When you are looking for a compact power wheel chair that is fun to drive and easy to transport – look no further than the Golden LiteRider Envy!



<b>Model Number</b>	<b>GP162</b>	
Medicare Code	Not Coded	
Weight Capacity	300 lbs.	
Drive Wheels	Rear	
Maximum Speed	3.5 mph	
Operating Range <sup>1</sup>	15.5 miles	
Ground Clearance	2.5"	
Turning Radius	28.5"	
Type Batteries	2-22 AH	
Colors	Cherry Red, Paradise Blue, Solar Flare Orange, Envy Green, Satin Silver, Sunburst Yellow	
Freewheel Mode	Yes	
Electronic Speed Control	Yes	
Electro-Mechanical Brakes	Yes	
Charger	Off Board	
Controller	40 Amp Dynamic Linx	
<b>Dimensions:</b>		
Length	36"	
Width	22.5"	
Height (ground to top of back of seat)	33" - 34"	
Height (ground to top of back of headrest)	N/A	
Ground to Top of Seat	20.5" - 21.5"	
Deck to Top of Seat	N/A	
Number of Seat Height Positions	2	
Size of Increments	1"	
<b>Weight of Unit: (assembled)</b>		
Front Half	37 lbs.	
Rear Half	36 lbs.	
Seat w/Arms	23 lbs.	
Batteries	27 lbs/pack	
<b>Tires:</b>		
Drive Wheels	9"	
Casters	6"	
Anti-Tip Wheels	3"	
<b>Standard Seat:</b>		
Stadium Style		
Back Height - no headrest	16"	
Back Height with headrest	N/A	
Width x Depth	17" x 17"	
Color	Black Vinyl	
<b>Optional Seats:</b>		
Stadium Style	Premium	
Back Height - no headrest	16"	17"
Back Height with headrest	N/A	N/A
Width x Depth	20" x 18"	17" x 16"
Color	Black Vinyl	Black Vinyl
<b>Seat Options:</b>		
Pneumatic Rise Seat Option	N/A	
Electric Rise Seat Option	N/A	
<b>Footrest:</b>		
Yes		
Height adjustable	Yes	
Angle adjustable	Yes	
<b>Warranty:</b>		
Frame	5 Years	
Drive Train	13 Months	
Electronics <sup>2</sup>	13 Months	
Batteries	12 Months	

<sup>1</sup> Battery range tested at 200 lbs. and will vary due to rider weight, drive surface, terrain and battery type.

<sup>2</sup> Electronics warranty excludes batteries.



Watch our video to learn about the benefits of the LiNX controller!

Scan here with smartphone camera.



Disassembles easily for transport!



Ergonomic LiNX controller



**NEW!** This model shown with Premium Seat Option!

### Available Accessories:

- Cane Holder Tube
- Oxygen Tank Holder
- Cup Holder
- Quad Cane Holder
- Walker Holder
- Crutch Holder
- Safety Flag
- Power Chair Cover
- Pack N' Go
- Seat Belt (Optional)



Perfect for transport!

Some accessories require a mounting bracket or mounting clips.