



Seat Color
Black Vinyl



Body Colors
Blue
Red



3 wheel GB117H
4 wheel GB147H

The Buzzaround XL HD features compact portability with a heavy-duty weight capacity! The Buzzaround XL HD features a robust 325 lb. weight capacity yet disassembles quickly and easily for transport. Standard features include a new LED angle-adjustable headlight, LED battery gauge, and more foot room than other scooters in its class. It's tight turning radius makes it ideal for indoor use, while it's ample ground clearance makes it great for outdoor use, as well.





LIGHTWEIGHT PORTABLE SCOOTER

Wireless disassembly in 4 easy steps!

| Item | GB117H | GB147H |
|--|-------------------------|-------------------------|
| Description | XLHD | XLHD |
| Medicare Code | K0801 | K0801 |
| Weight Capacity | 325 lbs. | 325 lbs. |
| Drive Wheels | Rear | Rear |
| Maximum Speed | 4 mph | 4 mph |
| Operating Range up to ¹ | 14 miles | 14 miles |
| Under Deck Clearance | 3" | 3" |
| Turning Radius | 37" | 47" |
| Type Batteries | 2-22 AH | 2-22 AH |
| Colors | Red, Blue | Red, Blue |
| Freewheel Mode | Yes | Yes |
| Electronic Speed Control | Yes | Yes |
| Electro-Mechanical Brakes | Yes | Yes |
| Charger | Off-board | Off-board |
| Controller | 50 Amp Dynamic R-Series | 50 Amp Dynamic R-Series |
| Dimensions | | |
| Length | 42.5" | 41.5" |
| Width | 21.5" | 21.5" |
| Height (ground to top of back of seat) | 33"-35" | 33"-35" |
| Ground to Top of Seat | 21"-23" | 21"-23" |
| Ground to Top of Deck | 4.75" | 4.75" |
| Deck to Top of Seat | 16"-18" | 16"-18" |
| Number of Seat Height Positions | 3 | 3 |
| Size of Increments | 1" | 1" |
| Front Axle to Rear Axle | 33.5" | 32" |
| Weight of Unit: (assembled)² | | |
| Front Half | 111 lbs. | 119 lbs. |
| Rear Half | 35 lbs. | 43 lbs. |
| Seat w/Arms | 25 lbs. | 25 lbs. |
| Batteries | 22 lbs. pack | 29 lbs. pack |
| Tires: | | |
| Front | Solid | Solid |
| Rear | 8" | 8" |
| Rear Anti-Tip Wheels | 2" | 2" |
| Standard Seat | | |
| Seat Back Height - no headrest | Stadium Style | Stadium Style |
| Width x Depth | 16" | 16" |
| Color | 18" x 16" | 18" x 16" |
| | Black Vinyl | Black Vinyl |
| Optional Seat: | | |
| Back Height - no headrest | Stadium Style | Stadium Style |
| Width x Depth | 16" | 16" |
| Color | 20" x 18" | 20" x 18" |
| | Black Vinyl | Black Vinyl |
| Warranty: | | |
| Frame | Lifetime Warranty | Lifetime Warranty |
| Drive Train | 2 Year Warranty | 2 Year Warranty |
| Electronics ² | 2 Year Warranty | 2 Year Warranty |
| Batteries | 1 Year Warranty | 1 Year Warranty |

¹. Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain, and battery type.

². Electronics warranty excludes batteries.



Step 1: Pop off the seat!



Step 2: Pull out the battery pack!



Step 3: Fold down tiller and lock in place!



Step 4: Pull forward on the drivetrain release lever to separate the front and rear frame sections!



Disassembles into 5 pieces!

4 wheel model shown



A battery charging port is conveniently located on the tiller under the control panel.

Features

- 325 lb. weight capacity
- Tiller mounted basket
- Angle-adjustable LED headlight
- 22 amp battery