



THE MOST STABLE BASE IN ITS CLASS!

YOU CAN'T STAND LIKE THIS  
ON ANY OTHER POWER WHEELCHAIR BASE

The Golden Compass™ line of power wheelchairs provide superior forward stability making them the safest power wheelchairs to mount or dismount. The Compass line has fewer moving parts than a power base with springs for the best forward stability and no tipping!

THE INNOVATIVE SEAT  
DESIGN MAKES OUR  
POWER WHEELCHAIRS  
THE EASIEST TO TRANSPORT!

Golden's unique Captain's Seat design found on the Compass Sport and Compass HD features the patented Fold-Flat Seat Bracket standard. This unique bracket allows the backrest to fold completely flat against the seat, requiring a lower clearance height to load the unit into a vehicle. The Compass HD also features the ergonomic Slide Seat Bracket so the seat can be placed onto or removed from the seat frame for easy assembly or disassembly.



POWER WHEELCHAIR PRODUCT SPECIFICATIONS

Model Number	GP605 Compass Sport	GP620 Compass HD	GP162 LiteRider Envy	GP205 Alanté Sport	GP207 Alanté
Medicare Code	K0823/K0822	K0825/K0824	Not Coded	K0823/K0822	K0816/K0815
Weight Capacity	300 lbs.	450 lbs.	300 lbs.	300 lbs.	300 lbs.
Drive Wheels	Center	Center	Rear	Front	Front
Maximum Speed	4 mph	4 mph	3.5 mph	3.5 mph	3.5 mph
Operating Range <sup>1</sup>	19.5 miles	12 Miles	15.5 miles	16.5 miles	16.5 miles
Ground Clearance <sup>3</sup>	3.25"	3.5"	2.5"	1.75"	1.75"
Turning Radius	19.5"	20"	26"	25"	25"
Type Batteries	2-U1	2-NF22	2-22 AH	2-U1	2-U1
Colors	Red, Blue	Red, Blue	Red, Blue	Red	Red
Freewheel Mode	Yes	Yes	Yes	Yes	Yes
Electronic Speed Control	Yes	Yes	Yes	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes	Yes	Yes	Yes
Charger	Off Board	Off Board	Off Board	Off Board	Off Board
Controller	Dynamic Linx 40 Amp	Dynamic Shark 75 Amp Remote	PG VR-2 50 Amp Remote	PG VR-2 50 Amp Remote	PG VR-2 50 Amp Remote
Length	39"	40"	33"	39.5"	39.5"
Width	24"	25"	22.5"	22.25"	22.25"
WARRANTY	Frame	Lifetime	5 Years	5 years	5 years
	Drive Train	13 months	13 months	13 months	13 months
	Electronics <sup>2</sup>	13 months	13 months	13 months	13 months

1. Operating range will vary due to rider weight, drive surface and terrain.  
2. Electronics warranty excludes batteries.  
3. Measured at the lowest point mid-frame.

Note: Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



# Golden Power Wheelchairs

HIGHLY DEPENDABLE WITH MAXIMUM MANEUVERABILITY,  
SUPERIOR STABILITY AND GREATEST FUNCTIONALITY



GOLDEN

BR-POWER.101315



# Live the Golden Life

Test drive the industry's most comprehensive line of consumer power wheelchairs today!

## Experience Constant Contact.

Inclines and thresholds are no match for the Golden Compass Sport's true center wheel drive design with the patented springless articulating chassis!

This unique suspension system keeps all six wheels on the ground simultaneously, eliminating the rocking sensation when traversing inclines, giving all Compass™ models unparalleled curb and threshold climbing ability. The Compass's suspension system pushes the center drive wheel down when the front and rear casters are raised and pushes the front and rear casters down when the drive wheels are raised. No other power base design can match the constant wheel to ground contact of any Golden Compass power base!



### LiteRider<sup>®</sup> GOLDEN ENVY

- Model GP162
- **Retail Cash Power Chair!**
- Wt. Cap.: 300 lbs.
- Turning Radius: 26"
- Batteries: 2-22AH
- Single Post Seat Mount
- DISASSEMBLES FOR TRANSPORT
- Operating Range: 15.5 mi.



### GOLDEN Alanté SPORT

- Model GP205
- Medicare Codes: K0823 & K0822
- Wt. Cap.: 300 lbs.
- Turning Radius: 25"
- Batteries: 2-U1
- Single Post Seat Mount



### Alanté BY GOLDEN

- Model GP207
- Medicare Codes: K0816 and K0815
- Wt. Cap.: 300 lbs.
- Turning Radius: 25"
- Batteries: 2 U-1
- Single Post Seat Mount



### golden COMPASS<sup>™</sup> SPORT

- Model GP605
- Medicare Codes: K0823 & K0822
- Wt. Cap.: 300 lbs.
- Turning Radius: 19.5"
- Batteries U1
- Single Post Seat Mount
- Dynamic Linx Controller Standard



### golden COMPASS<sup>™</sup> HD

- Model GP620
- Medicare Codes: K0825 & K0824
- Wt. Cap.: 450 lbs.
- Turning Radius: 20"
- Batteries NF22
- Four Post Seat Mount with Slide Seat Bracket