

**New  
Power Cloud!**



PR512MLA  
Medium/Large



PR510SME  
Small/Medium



PR-510MLA  
Medium/Large



**MaxiComfort  
Cloud**  
PR510 & PR512

The Cloud is the world's most comfortable lift chair recliner which comes in two sizes, Small/Medium and Medium/Large. Featuring our custom bucket seat, this chair is the ultimate in relaxation and rejuvenation. The newest edition to the Cloud family, the PowerCloud, also features our exclusive PowerPillow articulating headrest and a modern back design.

**PR510-SME  
FABRICS**



-  Hazelnut
-  Sterling
-  Copper
-  Calypso
-  Shiraz
-  Coffee Bean\*

**PR510-MLA  
FABRICS**

-  Hazelnut
-  Sterling
-  Copper
-  Calypso
-  Shiraz
-  Buckskin\*
-  Coffee Bean\*

**PR512-MLA  
FABRICS**

-  Hazelnut
-  Sterling
-  Copper
-  Calypso
-  Shiraz
-  Coffee Bean\*

-  Coffee Bean\*
-  Bride\*

\* upcharge applies



# MaxiComfort Cloud

M A X I C O M F O R T ™ S E R I E S

	Cloud	Cloud	Power Cloud
<b>Model Number</b>	PR510-SME	PR510-MLA	PR512-MLA
<b>Size</b>	Small/Med	Med/Large	Med/Large
<b>Weight Capacity</b>	375 LBS	375 LBS	375 LBS
<b>Positioning</b>	MaxiComfort	MaxiComfort	MaxiComfort
<b>Chaise Seat</b>	YES	YES	YES
<b>Pocketed Coil Seat</b>	NO	NO	NO
<b>Back Type</b>	Biscuit	Biscuit	Seam
<b>Removable Back Option</b>	NO	NO	NO
<b>Head Rest Cover</b>	YES	YES	NO
<b>Arm Rest Covers</b>	NO	NO	NO
<b>Battery Back Up</b>	YES	YES	YES
<b>Overall Width</b>	33.5"	38"	38"
<b>Overall Height</b>	41"	45"	46"
<b>Width Between Arms</b>	18.5"	20"	19"
<b>Floor to Top of Seat</b>	19"	20"	20"
<b>Seat Depth</b>	19.5"	21.5"	21"
<b>Seat to Top of Back</b>	27"	29"	30"
<b>Seat to Top of Arm</b>	7.5"	7"	7"
<b>Seat Angle Raised</b>	15°	15°	15°
<b>Floor to Top of Foot Rest Reclined</b>	23"	24.5"	24.5"
<b>Distance Required from Wall Reclined</b>	23"	24"	28.5"
<b>Extension reclined from Back of Seat to End of Foot Rest</b>	43"	44"	43"
<b>Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism</b>	YES	YES	YES
<b>Warranty on Electrical Parts</b>	3 Yrs.	3 Yrs.	3 Yrs.
<b>Pro Rated Years 4 Thru 7</b>	YES	YES	YES
<b>Warranty on Mechanical Labor</b>	3 Yrs.	3 Yrs.	3 Yrs.
<b>Warranty on Electrical Labor</b>	1 Yr.	1 Yr.	1 Yr.
<b>Smart Tek System Standard</b>	YES	YES	YES
<b>Customizable</b>	NO	NO	NO
<b>Left Hand Control Option</b>	YES	YES	YES
<b>Extra Pocket Option</b>	Standard	Standard	Standard
<b>Heat &amp; Massage Option</b>	YES	YES	YES
<b>Power Pillow Option</b>	NO	NO	Standard
<b>Power Adjustable Lumbar Option</b>	NO	NO	NO
<b>AutoFlex Option</b>	NO	YES	NO
<b>Foot Rest Extension Option</b>	NO	YES	NO
<b>Straight Lift Option</b>	NO	NO	NO
<b>USB Charging Port</b>	NO	NO	NO
<b>Place of Origin</b>	USA	USA	USA



## Golden Craftsmanship

Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.



## Added Comfort and Safety

Designed to perform year after year using advanced technology backed by the best warranty.



## Smart Engineering

Smooth, quiet lifting system offers a secure lift for the best performance.



## Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.



## Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.

\* See Customize Chart for Dimension limits  
 Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with  
 UL 73 and  
 CSA C22.2 No. 68

