

H E A V Y D U T Y P O V



Seat Color
Charcoal Vinyl

Body Colors
Orange
Red

Ride

Golden
Patriot™
4 wheel

If you're looking for the ultimate outdoor driving experience, look no further than the Golden Patriot™. The Golden Patriot's rugged outdoor design features an off road package, four-wheel independent suspension, a full lighting package and a full sized captain's seat. With a 400 lb. weight capacity and 13" solid, low profile tires, the Golden Patriot can give you the ride of your dreams.



Rise. Roll. Ride.

GOLDEN
technologies

YOUR LIFE IN MOTION

GR-575

Patriot 4 wheel Specifications

| | |
|--|--------------------|
| Model | GR575 |
| Medicare Code | Not Coded |
| Weight Capacity | 400 lbs |
| Drive Wheels | Rear |
| Maximum Speed | 8 mph |
| Operating Range ¹ | 13 miles |
| at 200 lbs | 18 miles |
| Ground Clearance | 3.5" |
| Turning Radius | 76.5" |
| Type Batteries | 2-75 AH |
| Colors | Red, Orange |
| Freewheel Mode | Yes |
| Electronic Speed Control | Yes |
| Electro-Mechanical Brakes | Yes |
| Mechanical Hand Brake | Yes |
| Charger | Off Board |
| Controller | PG S-140A |
| Dimensions: | |
| Length | 57" |
| Width | 28" |
| Height (ground to top of back of seat) | 46.5" - 48.5" |
| Height (ground to top of back of headrest) | 54" - 59" |
| Ground to Top of Seat | 27" - 29" |
| Ground to Top of Deck | 8.75" |
| Deck to Top of Seat | 18" - 20" |
| Number of Seat Height Positions | 3 |
| Size of Increments | 1" |
| Front Axle to Rear Axle | 40" |
| Weight of Unit: (assembled) | 314 lbs |
| Base | 160 lbs |
| Seat w/Arms | 50 lbs |
| Batteries | 52 lbs ea |
| Tires: | Solid, low profile |
| Front | 13" |
| Rear | 13" |
| Standard Seat: | Captain's Highback |
| Slide Seat | N/A |
| Back Height - no headrest | 21" |
| Back Height with headrest | 28.5" - 31.5" |
| Width x Depth (adjustable depth is standard, set to 18") | 20" x 18"/20" |
| Color | Charcoal Gray |
| Warranty: | |
| Frame | 5 Years |
| Drive Train | 1 Year |
| Electronics ² | 1 Year |

¹ Battery range will vary due to rider weight, drive surface and terrain.

² Electronics warranty excludes batteries.



Full front and rear lighting package standard



Floor deck features rugged rubber mat and anti-skid material



Four-wheel independent suspension standard



Water bottle holders for cold beverage bottles only. Removable for easy cleaning.



Control panel features safety horn, flashers, lights, turn signals, speed knob and backlit battery gauge



Black, low profile, non-marking solid tires with chrome hubcaps standard



Hand brake standard