

T R A D I T I O N A L S E R I E S



Rise

-  **W** WAREHOUSE
Hazelnut
-  **W** WAREHOUSE
Sterling
-  **W** WAREHOUSE
Copper
-  **B** BUILDERS
Calypso
-  **B** BUILDERS
Shiraz
-  **B** BUILDERS
Pearl

Cambridge
PR-401-MLA

Cambridge
PR-401-SME

Cambridge The popular Cambridge lift chair has been improved! The Cambridge now comes with Golden's exclusive hybrid sizing for a better fit. The new Small/Medium Cambridge fits people 5'1" to 5'6" and the new Medium/Large fits people 5'7" to 6'4". The new Cambridge features a more supportive biscuit backrest design and a softer, more comfortable seat.



Rise. Roll. Ride.

YOUR LIFE IN MOTION



Cambridge



	Cambridge PR-401-SME	Cambridge PR-401-MLA
Model Number	PR-401-SME	PR-401-MLA
Size	Small/Med	Med/Large
Weight Capacity	375 Lb.	375 Lb.
Recline Positions	3	3
Standard Fabrics	6	6
Fabric Options	All	All
Chaise Seat	YES	YES
Coil Springs	NO	NO
Back Type	Biscuit	Biscuit
Head Rest Cover	YES	YES
Arm Covers	YES	YES
Battery Back Up	YES	YES
Overall Width	33"	35"
Overall Height	41.5"	44.5"
Width Between Arms	20"	22"
Floor to Top of Seat	19"	22"
Seat Depth	19"	22"
Seat to Top of Back	27"	31"
Seat to Top of Arm	6.5"	7"
Seat Angle Raised	20°	20°
Floor to Top of Foot Rest Reclined	22"	23"
Floor to Top of Back Reclined	33"	35"
Back Angle Closed	105°	105°
Back Angle Reclined	140°	140°
Foot Rest Angle Closed	90°	90°
Foot Rest Angle Opened	180°	180°
Distance Required from Wall Reclined	18"	19"
Extension reclined from Back of Seat to end of Foot Rest	41.5"	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES
Warranty on Electrical Parts	3 Yrs.	3 Yrs.
Pro Rated Years 4 Thru 7	YES	YES
Warranty on Mechanical Labor	3 Yrs.	3 Yrs.
Warranty on Electrical Labor	1 Yr.	1 Yr.
Removable Back Standard	YES	YES
Smart Tek System Standard	YES	YES
Customizable*	NO	NO
Left Hand Control Option Available	YES	YES
Extra Pocket Option Available	YES	YES
Heat & Massage Option Available	YES	YES
Foot Rest Extension Option Available	YES	YES
Straight Lift Option Available	YES	YES
User Height (suggested)	5'1" - 5'6"	5'7" - 6'4"



501 M shown

Golden Craftsmanship

Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.



501 M shown

Added Comfort and Safety

Designed to perform year after year using advanced technology backed by the best warranty.



501 M shown

Smart Engineering

Smooth, quiet lifting system offers a secure lift for the best performance.



Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.



Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.

* See Customize Chart for Dimension limits
 Sizes may vary depending on fabric, filling material, upholstery and carpet thickness.
 All measurements taken on concrete floor using levels and metal rulers.