

Neb-u-Tyke™ Neb-a-Doodle

JOHN BUNN
A Graham-Field® Family Brand

This interactive, child-friendly design provides an advanced solution for effective aerosol treatment for children. Designed with an integrated drawing board, the Neb-a-Doodle will keep child relaxed and entertained during aerosol treatment and can be used over and over again!

The John Bunn Neb-u-Tyke Neb-a-Doodle is the latest addition to the Graham-Field line of quality nebulizer compressors that offer innovation, performance, and value for both patients and providers.



Features:

- Reliable, piston-driven compressor
- Interactive, drawing board for child to draw, doodle and write during aerosol treatment
- Convenient and compact design for easy transport and storage
- Comes complete with Nebulizer Accessory Kit
- Five year limited warranty
- Designed to meet the requirements of HCPCS Code: E0570

JB0112-070 1 ea

graham
field

Neb-u-Tyke™ Neb-a-Doodle



- The John Bunn Neb-a-Doodle Nebulizer comes complete with:
- 1 Neb-u-Tyke Neb-a-Doodle Nebulizer Compressor
 - 1 User Manual
 - 1 Nebulizer Accessory Kit complete with:
 - 1 Reusable nebulizer (single patient use only)
 - 1 Adult mask
 - 1 Pediatric mask
 - 1 Mouthpiece
 - 1 PVC tubing
 - 1 "T" adapter
 - 3 Spare filters
 - 1 Reservoir tube

Technical Specifications

Power	115VAC, 60Hz
Power Consumption	≤ 60 watts
Operating Temperature Range	41°F TO 104°F
Operating Humidity Range	30 to 95% RH
Storage Temperature Range	-40°F TO 158°F
Storage Humidity Range	10 to 95% RH
Dimensions	10" x 6.5" x 4.2" (255 mm x 166mm x 108 mm)
Weight	4 lb
Max Pressure (Compressor)	32.5 psi
Max Flow (Compressor)	Approx 12 LPM
Particle Size (Nebulizer)	≤5 microns
Capacity (Nebulizer)	10 ml
Method Of Nebulization	Piston Driven
Sound Level	under 57 dBA

graham
field