

Neb-u-Tyke Penguin Nebulizer Compressor

- Small footprint and light weight provide exceptional portability
- Compressor weighs only 3.5 lbs
- Reliable, Piston-driven Compressor
- Colorful, Child-Friendly
- 5 year limited warranty
- Designed to meet the requirements of HCPCS Code: E0570

JB0112-062	Penguin Nebulizer	1 ea
JB0112-06F	Filter for Penguin Nebulizer	1 ea
JB77100	Penguin Nebulizer And Igloo	1 ea
JB76210	Igloo case only	1 ea
IGLOO-DISPLAY	Retail display	1 ea

The Neb-u-Tyke Penguin Nebulizer comes complete with:

- 1 disposable nebulizer (single patient use only)
- 1 mouthpiece
- 1 pediatric mask
- 1 "T" adapter
- 1 PVC tubing
- 3 spare filters
- 1 reservoir tube



retail display



with optional igloo carrying case



Technical Specifications

Method Of Nebulization	Piston Driven
Power	115VAC, 60Hz
Power Consumption	63 watts
Operating Temperature Range	41°F to 104°F (5°C to 40°C)
Operation:	30% to 95% RH
Storage Temperature Range	-40°F to 158°F (-40°C to 70°C)
Storage:	30% to 95% RH

Weight	3.5 lb (1.59 kg)
W x D x H	6.7" x 4.8" x 6.5" (17.0 cm x 12.2 cm x 16.5 cm)
Max Pressure (Compressor)	32.5 psi
Flow Range (Compressor)	0 - 10 liters/min
Particle Size (Nebulizer)	0.5 - 5 microns
Capacity (Nebulizer)	10 ml