


The Everest & Jennings Bariatric Transport Chair is designed as an affordable patient transfer device. Increased weight capacity and a wider width offer the user increased strength, comfort and durability.



Features:

- Fold down back makes transport and storage easy
- Black, padded, Nylon upholstery is flame retardant
- 24" seat width
- Carry straps on seat allow for ease of lifting and folding
- Desk length padded arms provide user comfort
- Weighs 31 lb (w/o front rigging)
- Maximum Weight Capacity 400 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of Medicare code E1039

Warranty:

Three (3) year limited warranty



Bariatric Transport Chair		Dimensions and Weights								
		Seat			Back	Overall			Wt.	
EJ777-3	Aluminum - Silver	24"	18"	20"	18"	40"	39 1/4"	30 1/4"	12"	31 lb

All dimensions are ±1/4"
*Without Front Rigging