



SUPPORT SURFACES

BARIATRIC SKIN PROTECTION CUSHIONS

The Bariatric Skin Protection Cushions provides pressure distribution and weight equalization to meet unique demands of the bariatric patient.



Cushion Features

- 4" cushion
- Medium density foam base
- Non-aqueous specialized visco-elastic gel core layer reduces tissue interface pressure at critical points
- Proprietary gel is a flexible solid that will not leak
- Gel is placed at three pressure zones
- Visco-elastic memory foam topper for superior comfort and pressure relief
- Positions the user with a moderate slope from front to back

Cover Features

- Fluid-resistant stretch top cover, reduces shear and friction
- Antimicrobial
- Flame retardant
- Removable and washable
- Non-skid base
- Strap to secure cushion firmly to the chair



Specifications

Item	UOM	Meets HCPCS	Overall Width		Overall Length		Overall Thickness		Weight Capacity EVENLY DISTRIBUTED		Made in USA	Warranty in months
			inch	cm	inch	cm	inch	cm	lb	kg		
20252218B	1/ea	E2604	22	55.9	18	45.7	4	10.2	650	294.8	■	18
20252418B	1/ea	E2604	24	61	18	45.7	4	10.2	650	294.8	■	18
20252619B	1/ea	E2604	26	66	19	48.3	4	10.2	650	294.8	■	18
20253019B	1/ea	E2604	30	76.2	19	48.3	4	10.2	650	294.8	■	18



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.



Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

© 2008, GF Health Products, Inc. All Rights Reserved.
Gendron is a registered trademark of GF Health Products, Inc.