

## Labstar™ Deluxe Sphygmomanometer

- Accumax™ nylon cuff with index & range markings
- Latex two (2) tube bladder
- Large knurled air release valve for precision deflation
- Deluxe end valve reduces dust build-up
- Color coordinated gauge and cuff - blue, black, grey or orange
- Pressure measurement from 20 mmHg to 300 mmHg
- Lifetime calibration warranty
- Zippered carrying bag is included

**Caution:** This product contains natural rubber latex, which may cause allergic reactions.

<b>202</b>	Adult, Blue	1 ea
<b>202C</b>	Child, Blue	1 ea
<b>202I</b>	Infant, Blue	1 ea
<b>202N</b>	Newborn, Blue	1 ea
<b>202T</b>	Thigh, Blue	1 ea
<b>202X</b>	Large Adult, Blue	1 ea
<b>202BK</b>	Adult, Black	1 ea
<b>202CBK</b>	Child, Black	1 ea
<b>202IBK</b>	Infant, Black	1 ea
<b>202TBK</b>	Thigh, Black	1 ea
<b>202XBK</b>	Large Adult, Black	1 ea
<b>202GY</b>	Adult, Grey	1 ea
<b>202XGY</b>	Large Adult, Grey	1 ea
<b>202COR</b>	Child, Orange	1 ea
<b>202IOR</b>	Infant, Orange	1 ea
<b>202OR</b>	Adult, Orange	1 ea
<b>202XOR</b>	Large Adult, Orange	1 ea



**202X**



**202GY**



**202BK**



**202**



**202T**

<b>Cuff Class</b>	<b>Limb Range (inches)</b>	<b>Bladder (inches)</b>	<b>Cuff Dimensions (inches)*</b>
Newborn	3.75" - 5.25"	1" x 3.55"	7" W x 2" H
Infant	5" - 7.5"	2.35" x 4.72"	10.5" W x 3" H
Child	7" - 10.5"	3.55" x 5.90"	13.5" W x 4" H
Adult	9.5" - 15"	4.72" x 8.66"	20.5" W x 5.5" H
LG Adult	13" - 19"	5.71" x 12.41"	25" W x 6.5" H
Thigh	17" - 26"	7.09" x 14.57"	29.5" W x 8.5" H

\* These are subjective measurements and can vary ± 0.25"

