

Bariatric Transfer Bench

- Provides assistance to those who have difficulty stepping over bathtub walls
- Durable blow-molded plastic bench and backrest provide a 600 lb. weight capacity
- Drainage holes in seat minimize water build-up
- Extra-deep seat provides added stability
- Anodized aluminum frame is lightweight, durable and rust-resistant
- Seat height is adjustable in 1" increments
- Reversible back attachment to accommodate any bathroom
- Attractive platinum grey color
- Comes assembled with tool-free back attachment
- Maximum Weight Capacity 600 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: EO248
- Lifetime limited warranty

7925A 1 ea



Specifications

Overall Width	29 1/2"
Overall Depth.....	20 1/2"
Seat Width	26 3/4"
Seat Depth.....	16 1/4"
Minimum Seat Height	16"
Maximum Seat Height	20"
Clearance Beneath Frame	14" - 18"
Height of Back Support Above Seat	13"
Maximum Weight Capacity.....	600 lb EVENLY DISTRIBUTED

