

## Deluxe Shower Transport Chair with Footrests

- Designed for safe transport to and from the shower
- Anodized, non-corrosive aluminum frame construction will not rust
- Comfortable nylon backrest removes for easy cleaning
- Cushioned armrests enhance user comfort
- 5" swivel casters with rear wheel safety locks facilitate maneuverability in tight spaces
- Rolls easily over most toilets
- Allows complete access to body for easy cleaning during bathing
- Powder-coated steel front brace with footrest
- Foot plates tip up and out of the way for unobstructed entry and exit
- Footrests comfortably support the feet during transport
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Limited lifetime warranty

7915A-1      1 ea

### Specifications

Width Between Arms .....	17 1/2"
Overall Width .....	21 1/2"
Overall Depth.....	30 1/2"
Overall Height.....	38"
Seat Width.....	14"
Seat Depth.....	16 1/2"
Seat Height.....	22"
Maximum Weight Capacity .....	250 lb EVENLY DISTRIBUTED

