

Platinum Collection Deluxe Padded Commode Bath Seats

- Open front seat provides easier access for personal care
- Comfortable, cushioned seat and backrest are water-tight and easy to clean
- Anodized aluminum frame is lightweight, durable and rust-resistant
- Seat height is adjustable in 1" increments
- Attractive platinum grey color
- Comes complete with commode pail and cover
- Available with or without support arms
- Comes assembled with tool-free arm and back attachments
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0245
- Lifetime limited warranty

7946A without support arms 1 ea

7947A with support arms 1 ea

<u>Specifications</u>	<u>7946A</u>	<u>7947A</u>
Overall Width.....	20 ½"	20 ½"
Overall Depth.....	21"	21"
Seat Width.....	15 ¾"	15 ¾"
Seat Depth.....	16"	16"
Minimum Seat Height.....	17 ½"	17 ½"
Maximum Seat Height.....	22 ½"	22 ½"
Height of Back Support Above Seat.....	16"	16"
Diameter Commode Opening.....	8"	8"
Width Between Arms.....	NA	17 ¾"
Maximum Weight Capacity.....	300 lb EVENLY DISTRIBUTED	

