

# When **Bigger IS** Better...

## Introducing the new Lumex Set N' Go Wide

The Set n' Go Wide Height Adjustable Rollator offers the convenience of a two-in-one product all rolled into one with added features for larger individuals. With its enhanced, wider spacing

of 19" in between the handgrips and a larger seat, the Set n' Go Wide is ideal for those users needing a little more room for comfortable movement and seating.



19" opening width between the handgrips is designed for larger individuals ...no more squeezing into a smaller product



**WIDE** 15.75" seamless seat provides extra padding and seating area for added comfort to support 350 lb weight capacity **EVENLY DISTRIBUTED**

Inner leg tubing is precisely engineered to reduce instability and secured with leg knob for added reinforcement



Large 8" casters for indoor and outdoor terrain



Velcro folding strap keeps rollator folded for easy storage and transport



Easy height adjustment with clear position indicators to accommodate a broad range of user heights

Available in Blue (RJ4718B), Burgundy (RJ4718R) and Silver (RJ4718S)

**Specifications\***


Width Inside Handgrips:	19.4" (492 mm)
Overall Depth:	26.9" (683 mm)
At Minimum Seat Height:	27.5" (699 mm)
At Maximum Seat Height:	31.5" (800 mm)
<b>Height</b>	
At Minimum Handle Height:	32" (813 mm)
At Maximum Handle Height:	38" (965 mm)
<b>Seat</b>	
Height:	19.5" – 23.5" (495 mm - 597 mm)
Width:	15.75" (400 mm)
Depth:	12" (305 mm)
Casters Diameter:	8" (203 mm)
Maximum User Weight:	350 lb (159 kg) Evenly Distributed
Product Weight:	16.5 lb (7.5 kg)

\*Specifications are subject to change without notice



**WIDEN YOUR OFFERING TO YOUR CUSTOMERS AND  
GET A SET N' GO WIDE TODAY!**

**More Lumex Offerings**



**GRAHAM-FIELD**