

Versa-Helper Trapeze

- Helps patient safely change position in bed and transfer from bed to other equipment with minimal assistance.
- Rubber bumper guards protect walls.
- Handbar location fully adjustable.
- Height adjustable from 43" to 52".
- VersaGuard-coated clamps and handbar.
- For use with metal frame headboards only.
- Assisted weight capacity of 250 lbs.
- Designed to meet the requirements of HCPCS code: E0910.

2800A	Chrome-plated	1 ea
2800GA	Gray	1 ea



2800GA



2800GA Trapeze Versa-Helper with
2840GA floor stand.

