

Lumineb II

The Lumineb II Nebulizer Compressor is an all included solution when you need to administer prescribed medication via nebulization

Features include:

- Lightweight, high-quality, low-noise nebulizer
- Compact design with rounded corners with integrated carrying handle for transport
- Piston-driven motor ensures high performance and long life
- Easy to clean and maintain
- Designed to meet the requirements of HCPCS code: E0570
- Five year limited warranty

Each nebulizer comes with:

- (1) Disposable nebulizer cup
- (1) Mouthpiece
- (1) 7' Clear Vinyl tubing
- (5) replacement filters

Optional Accessories Available:

- Adult Mask (Item 5704)
- Pediatric Mask (Item 5703)
- Carrying Case (Item 5703CC)

5710 1 ea

Specifications

Method of Nebulization:	Piston Driven
Power:	120V, 60 Hz
Operating Temperature:	50F to 104F (10C to 40C)
Operating Humidity Range:	10-95% RH
Storage Temperature Range:	-4F – 158F (-20C to 70C)
Storage Humidity Range:	10-95% RH
Weight:	3.3 lbs (1.5 kg)
Width:	7.75" (195mm)
Height:	6" (150mm)
Length:	8.5" (210mm)
Compressor Pressure Range:	30- 36 PSI (210 – 250 kPa)
Operating Pressure Range:	.8 – 16 PSI (50 – 100 kPa)
Flow Range:	6 -8 LPM
Particle Size:	4um
Average Nebulization Rate:	0.2ml/min
Capacity (nebulizer):	5ml
Sound Level:	60 dBA

