

Multi-position Open Padded Raised Toilet

- Unique open seat promotes patient independence and facilitates perineal cleaning and bowel training
- Front opening is 4 ¾" wide
- Seat can be rotated from front to either side or rear mounting for maximum access to perineal area
- Thick, heat-sealed padded seat is ideal for patients with circulatory problems or patients who need additional comfort
- Angle and height are adjustable to maximize comfort and reduce bending and patient effort required
- Comes complete with four heavy duty VersaGuard coated locking brackets
- Fits most toilets
- Optional extra-wide bracket (6499A) available for extra-wide toilet bowls
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0244
- Lifetime limited warranty

6497A 1 ea



Specifications

Seat Width	16"
Seat Depth	16"
Seat Height Range	6" - 7 ½"
Maximum Weight Capacity	300 lb EVENLY DISTRIBUTED



GRAHAM-FIELD