

The Paramount XD is ideal for the bariatric user. It features 26" and 30" seat widths as well as a weight capacity of 650 lbs. Meets Medicare specifications for K0007 for wheelchairs.



Standard Features

- Black reinforced nylon upholstery
- 19 1/2" Seat depth
- Swingaway footrest and Elevating legrest options
- Maximum weight capacity 650 lb - EVENLY DISTRIBUTED
- Reinforced welds and gussets on frame
- 8"x 2" Solid casters
- 24" Solid rear wheels

Warranty

Five (5) year limited warranty

Paramount XD		Dimensions and Weights								
		Seat			Back	Overall				Wt.
		Icon	Icon	Icon	Icon	Icon	Icon	Icon	Icon	Icon
5PX10620	Detachable Desk Arm, Swingaway Footrest	26"	19 1/2"	20 1/2"	18 1/2"	49"	39"	35 1/2"	14 1/2"	93 lb
5PX10630	Detachable Desk Arm, Elevating Legrest	26"	19 1/2"	20 1/2"	18 1/2"	49"	39"	35 1/2"	14 1/2"	93 lb
5PX10820	Detachable Desk Arm, Swingaway Footrest	30"	19 1/2"	20 1/2"	18 1/2"	49"	39"	40"	14 1/2"	93 lb
5PX10830	Detachable Desk Arm, Elevating Legrest	30"	19 1/2"	20 1/2"	18 1/2"	49"	39"	40"	14 1/2"	93 lb