

**Splash™ Bath Lift**

**An ultra-lightweight, budget-friendly bath lift**

With all the basic elements needed to safely maintain personal hygiene, the Splash bath lift has been designed to offer good value for the budget conscious. Weighing an incredibly low 14.3 lbs, Splash is not only the lightest battery-powered bath lift available in the US, its robust design also comfortably lifts users of up to 280 lbs.

**A simple, practical bathing solution**

**Fixed back design**

- Optimizes legroom by allowing Splash to fit close to the back of a tub

**Low seat travel**

- Maximizes water immersion for improved personal hygiene

**Smooth movement and solid, stable design**

- Offers a safe bathing experience

**Slim remote with large hand controls for easy operation**

- Can be easily operated by users with small hands and/or limited dexterity

**Battery level indicator with Stop Descent™ feature**

- Monitors battery power, preventing descent when low

**Simple to clean surfaces with limited crevices**

**Trim and lightweight for easy moving and handling**

**Total weight of only 14.3 lbs**

- Offers quick and uncomplicated fitting and removal from the bath in one complete unit

**Ultra-compact design**

- Provides ultimate ease of carrying
- Can be stored in even the smallest of spaces

**Subtle, unimposing, and unobtrusive design**

- Blends discreetly into any bathroom environment
- Barely noticeable when stored in its lowest position in the bath


**Slim and Modern look**

- Provides a fresh, non medical appearance in the home

5033A-1 1 ea

**Specifications**

Seat Width.....	13"	Seat Depth.....	13.4"
Flap Width.....	6.9"	Total Depth.....	20.1"
Total Width.....	26.8"	Net Weight (excluding battery).....	14.3 lbs
Total Height.....	22.8"	Shipping Weight.....	19.8 lbs



Lumex® is a registered trademark of GF Health Products, Inc. Splash™ is a trademark of Mountway Limited.  
Lumex Splash is powered by Mountway, LTD.

