


# Bathroom adaptations

LinidoSolutions

## Technical drawings / Door frame grab rails

All measurements are in mm. For complete article numbers, please consult the order information document which is available on our website.

**Door frame grab rail S/M/L** | Art. no. LI2619.0301-02, LI2619.0451-02, LI2619.0601-02



| Art. Nr    | A   | B   | C   | D   |
|------------|-----|-----|-----|-----|
| LI2619.030 | 300 | 280 | 240 | 215 |
| LI2619.045 | 450 | 430 | 390 | 365 |
| LI2619.060 | 600 | 580 | 540 | 515 |