

Standard Crib

Made in USA



HCPCS code E1399

Standard Features

- Manual Head & Knee Elevation
- 2 Drop Sides
- 3" Foam Mattress
- 5" Heavy Duty Casters (1 Steer, 2 Brake, 1 Swivel)
- Confetti Epoxy Finish
- Steel Construction

Options

- IV Pole- Add up to 4
- Undercrib Shelf
- Oxygen Tank Holder
- 5" Foam Mattress
- 5" Pressure Reduction Mattress
- Additional Epoxy Colors Available at No Additional Charge

Articulation

- Head to 35 degrees
- Thigh/Foot to 15 degrees

Our Standard Cribs are intended for use by a child who has limited mobility (cannot sit up, high kneel or stand) and will not be at risk of climbing over the rail and falling out.

HEMOCARE STANDARD CRIB				
MODEL #	H1942-CGP	H1952-KGP	H1973-MGP	H1982-YGP
Overall Dimensions	30" x 60"	36" x 60"	36" x 72"	36" x 83"
Mattress Size	28" x 58"	34" x 58"	34" x 70"	34" x 81"
Overall Crib Height	56"	56"	56"	56"
Floor to Top of Mattress Platform	31"	31"	31"	31"
Height of Side Rails	26"	26"	26"	26"

