




**2013** Outside Lift Systems





MADE IN USA



**AL100**  
Universal Scooter Lift

**AL160**  
Profile Scooter Lift







**AL300 / AL300HD**  
Fusion Lift

**AL300RV**  
RV Lift

**AL301XL / AL301XLHD**  
Fusion Lift

<b>Accommodates</b>	Scooters	Scooters	Power Chairs & Scooters	Power Chairs & Scooters	Larger Power Chairs & Scooters
<b>Platform Size</b>	27.5" x 39" aluminum deck, adjustable cradle extends up to 46"	27.5" x 23"	27.5" x 39" aluminum deck, adjustable cradle extends up to 46"	27.5" x 39" aluminum deck, adjustable cradle extends up to 46"	36" x 39" aluminum deck, adjustable cradle extends up to 53"
<b>Lifting Weight Capacity</b>	350 lbs. (160 kg)	350 lbs. (160 kg)	350 lbs. (160 kg) / 400 lbs.	350 lbs. (160 kg)	350 lbs. (160 kg) / 400 lbs.
<b>Installed Weights</b>	83 lbs. (38 kg)	74 lbs. (34 kg)	80 lbs. (36 kg) / 82 lbs. (37 kg)	82 lbs. (37 kg)	100 lbs. (45 kg) / 102 lbs. (46 kg)
<b>Maximum Wheel Base</b>	up to 42" wheelbase	up to 42"	up to 42" wheelbase	up to 42" wheelbase	up to 48" wheelbase
<b>Securement</b>	Automatic Hold-down Foot and Wheel Cradle	Automatic Hold-down Foot and Wheel Cradle	4 Q'Strain Retractors and Wheel Cradle	4 Q'Strain Retractors and Wheel Cradle	4 Q'Strain Retractors and 2 Wheel Cradles
<b>Hitch Class</b>	II or greater	I or greater	II or greater	III or greater	III or greater

<b>Benefits</b>	<ul style="list-style-type: none"> <li>One Switch Operation</li> <li>Quick, Automatic Securement</li> <li>No Scooter Modifications</li> <li>License Plate Mount Included</li> <li>Harmar Blue Touch Technology</li> <li>3-Year Warranty</li> </ul>	<ul style="list-style-type: none"> <li>Streamlined, Open Platform</li> <li>One Switch Operation</li> <li>Automatic Hold-down Foot</li> <li>Carries Scooter up to 350 lbs.</li> <li>No Scooter Modifications</li> <li>License Plate Mount Included</li> <li>3-Year Warranty</li> </ul>	<ul style="list-style-type: none"> <li>One Switch Operation</li> <li>Fits a Power Chair or Scooter</li> <li>No Scooter or Chair Modifications</li> <li>License Plate Mount Included</li> <li>Harmar Blue Touch Technology</li> <li>3-Year Warranty</li> </ul>	<ul style="list-style-type: none"> <li>One Switch Operation</li> <li>Engineered to Handle the Added Stress from RV's</li> <li>Heavy Duty Hitch Adapter Standard (26" Hitch Height) Upgrade Available</li> <li>No Scooter or Chair Modifications</li> <li>License Plate Mount Included</li> <li>3-Year Warranty</li> </ul>	<ul style="list-style-type: none"> <li>One Switch Operation</li> <li>Largest Platform in the Industry, a Harmar Exclusive</li> <li>No Scooter or Chair Modifications</li> <li>License Plate Mount Included</li> <li>Harmar Blue Touch Technology</li> <li>3-Year Warranty</li> </ul>
<b>Securement Methods</b>	 <p>Automatic Hold-down Foot</p>  <p>Automatic Hold-down Arm</p>  <p>Q'Strain Retractable Hooks</p>  <p>Dock 'N' Lock</p>				

					
<b>AL500 / AL500HD</b> Universal Power Chair Lift	<b>AL560</b> Automatic Universal Power Chair Lift	<b>AL560XL</b> Automatic Extra Large Power Chair Lift	<b>AL570 / AL570XL</b> Power Chair Lift with Dock 'N' Lock	<b>AL580</b> Profile Mid-Wheel Drive Lift	<b>AL580XL</b> Extra Large Profile Mid Wheel Drive Lift
Power Chairs	Power Chairs	Larger Power Chairs	Power Chairs / Larger Power Chairs	Mid-wheel Power Chairs with 10" Drive Wheels	Mid-wheel Power Chairs with 14" Drive Wheels
27.5" x 39" aluminum deck	28.5" x 38.5" adjustable platform	32.5" x 38.5" adjustable platform	26.5" x 34.5" / 35" x 34.5"	16.5" x 27"	17.5" x 27"
350 lbs. (160 kg) / 400 lbs.	350 lbs. (160 kg)	350 lbs. (160 kg)	350 lbs. (160 kg)	350 lbs. (160 kg)	350 lbs. (160 kg)
75 lbs. (34 kg) / 77 lbs. (35 kg)	89 lbs. (40 kg)	101 lbs. (46 kg)	69 lbs. (31 kg) / 79 lbs. (36 kg)	63 lbs. (29 kg)	63 lbs. (29 kg)
N/A	N/A	N/A		N/A	N/A
4 Q'Straint Retractors	Hold-down Arm	Hold-down Arm	Dock 'N' Lock	Hold-down Arm	Hold-down Arm
II or greater / III or greater	II or greater	III or greater	II or greater / III or greater	II or greater	II or greater
<ul style="list-style-type: none"> <li>One Switch Operation</li> <li>Patented, 4-strap Securement</li> <li>No Chair Modifications</li> <li>License Plate Mount Included</li> <li>Harmar Blue Touch Technology</li> <li>3-Year Warranty</li> </ul>	<ul style="list-style-type: none"> <li>One Switch Operation</li> <li>Adjustable Platform to Provide a Custom Fit</li> <li>No Chair Modifications</li> <li>License Plate Mount Included</li> <li>Harmar Blue Touch Technology</li> <li>3-Year Warranty</li> </ul>	<ul style="list-style-type: none"> <li>One Switch Operation</li> <li>Adjustable Platform to Provide a Custom Fit for Larger Power Chairs</li> <li>No Chair Modifications</li> <li>License Plate Mount Included</li> <li>Harmar Blue Touch Technology</li> <li>3-Year Warranty</li> </ul>	<ul style="list-style-type: none"> <li>Dock 'N' Lock automatically secures the power chair to the platform using a positive lock docking device.</li> <li>Adjusts to accommodate a wide range of power chairs.</li> <li>Swing-away option available</li> <li>License Plate Mount Included</li> <li>3-Year Warranty</li> </ul>	<ul style="list-style-type: none"> <li>One Switch Operation</li> <li>Compact and Lightweight</li> <li>No Chair Modifications</li> <li>License Plate Mount Included</li> <li>Harmar Blue Touch Technology</li> <li>3-Year Warranty</li> </ul>	<ul style="list-style-type: none"> <li>One Switch Operation</li> <li>Compact and Lightweight</li> <li>No Chair Modifications</li> <li>License Plate Mount Included</li> <li>Harmar Blue Touch Technology</li> <li>3-Year Warranty</li> </ul>

## Popular Accessories

WHEELCHAIR AND SCOOTER COVERS



HITCHES



SWING-AWAY



BATTERY PACK



## Also Available from Harmar

INSIDE LIFTS



HYBRID LIFTS



VERTICAL LIFTS



STAIR LIFTS



INCLINE LIFT



POOL LIFT



DUMBWAITER



HOME ELEVATORS

