

A. Model 1418
Bariatric Space-Saver
Therapy Platform



FEATURES:

- (6) High pressure laminate legs and a sturdy plywood sub-top.
- 2" High-density urethane foam top.
- 1000 lbs. patient weight capacity.
- Laminate: 346 Natural Oak (shown) or 927 Folkstone Gray.
- Upholstery: Choice of (9) Pro-Form™ Vinyl Colors.

MODEL	LENGTH	WIDTH	HEIGHT
1418	78"	30"	18"



A. Model 1418

B. Model 1422
Bariatric Wall Folding
Therapy Platform



FEATURES:

- Fold-away design saves floor space.
- 2" High-density urethane foam top.
- Sturdy, floor-supported wall mounting plate for added strength.
- Projects only 8" when folded.
- 1000 lbs. patient weight capacity.
- Laminate: 346 Natural Oak (shown) or 927 Folkstone Gray
- Upholstery: Choice of (9) Pro-Form™ Vinyl Colors.

MODEL	LENGTH	WIDTH	HEIGHT
1422	78"	30"	18"



B. Model 1422

C. S-201, S-202, S-205, S-206
Econo Mat Platforms



FEATURES:

- Sturdy panel end legs with wear-resistant Folkstone Gray laminate.
- 2" foam standard, 18" height.
- Upholstered in seamless nylon reinforced 113 Blue or 119 Gray.

MODEL	COLOR	WIDTH	LENGTH	HEIGHT
S-201	Blue	4'	7'	18"
S-202	Gray	4'	7'	18"
S-205	Blue	5'	7'	18"
S-206	Gray	5'	7'	18"



C. Model S-205