

## A. Model 4082 Combination Treatment/ Work Table/Desk



### FEATURES:

- Dual-purpose – it's a treatment table or a desk. Ideal for schools or clinics.
- 2" thick removable reversible urethane foam pad with nylon reinforced vinyl cover.
- Top pad has straps and can hang vertically on table when not in use.
- Folkstone Gray high pressure laminate sub-top with black 3mm edging.
- Long stretcher offset to the back for easy access when sitting.
- Custom laminate and upholstery color choices available for Contract Orders of 10 or more tables.

### COLORS:

**Laminate** – Folkstone Gray

**Upholstery** – Pad is Blue on one side and Gray on reverse side with 2" Gray border.

*(FOLD OUT FRONT COVER)*

MODEL	LENGTH	WIDTH	HEIGHT
4082	72"	24"	32"

## OPTIONS

FOR MODELS 4082, 4075

- \*00 – Special Table Height (for 4082 only)
- \*04 – Table Safety Strap
- \*07 – Nose Cutout (for 4075 only)
- 27 – Paper Dispenser & Cutter Combo
- 32 – Pillow, full size
- 35 – Pillow, small

\*FACTORY INSTALLED

SEE PAGE 4

## B. Model 4075 Wall Folding Treatment/ Changing Table



### FEATURES:

- Fold-away design saves floor space.
- 2" high-density foam upholstered top.
- Sturdy, floor-supported wall mounting plate for added strength.
- (2) folding high pressure laminate legs.
- Projects only 8" when folded.
- Wall fasteners not included.

### COLORS:

**Laminate** – Natural Oak

**Upholstery** – 17 standard vinyl colors

*(FOLD OUT FRONT COVER)*

MODEL	LENGTH	WIDTH	HEIGHT
4075	72"	24"	31"



**A. Model 4082  
"Pad On" — It's a  
Treatment Table**  
*Shown with optional  
#35 pillow*



**Model 4082  
"Pad Off" —  
It's a Desk**



**B. Model 4075**



**Model 4075 (folded)**