

Common Features

SPACE-SAVER TREATMENT TABLES —
MODELS 4102, 4124, 4143:



- Hausmann Space-Saver Tables meet your needs for a full size table in a compact space. When closed, tables are only 52" long. When open, they extend to a conventional 72".
- Full width 20" long foot end folds down to reduce table length to 52".
- Positive action locking device allows the foot end to be easily elevated into three positions: 55°, 75° and 90°.
- 2" High-Density urethane foam top.
- Resist-All™ vinyl upholstery finish.

COLORS:

Laminate – Natural Oak, Wild Cherry, Folkstone Gray

Upholstery – 17 standard vinyl colors
(FOLD OUT FRONT COVER)

OPTIONS

FOR 4100 SERIES TABLES

- *04 – Table Safety Strap
- *05 – Headrest, adjustable
- *06 – Backrest, adjustable
- *07 – Nose Cutout
- 23 – Paper Dispenser
- 25 – Paper Cutter
- 27 – Paper Dispenser & Cutter Combo
- 32 – Pillow, full size
- 35 – Pillow, small
- *99L – Left Hand Table Model

*FACTORY INSTALLED

SEE PAGE 4



A. Model 4143

A. Model 4143

SPECIAL FEATURES:

- Fully enclosed all laminate cabinet.
- (2) hinged doors with 3mm edging.
- Concealed, self-closing Euro-hinges.
- Laminate interiors with (1) adjustable shelf.
- Black laminated plywood base.



MODEL	LENGTH	WIDTH	HEIGHT
4143	52" - 72"*	24"	31"

* Fully extended



see Footstools
pages 33, 34

B. Model 4102

SPECIAL FEATURES:

- High pressure laminate legs, apron and stretchers with H-Brace.

MODEL	LENGTH	WIDTH	HEIGHT
4102	52" - 72"*	24"	31"

* Fully extended.

C. Model 4124

SPECIAL FEATURES:

- H-Brace with full width laminate shelf.
- High pressure laminate legs, apron and stretchers.

MODEL	LENGTH	WIDTH	HEIGHT
4124	52" - 72"*	24"	31"

* Fully extended.

B. Model 4102



C. Model 4124

Shown with option:
06 – Backrest

