

Common Features

MAXIMUM VALUE EXAM TABLES –
SERIES 4420, 4421, 4423, 4424:



- Removable vinyl top pad attached with VELCRO® brand fasteners.
 - Contoured padded head section.
 - Backrest is pneumatic Air-Spring with dual adjustment controls from 0° to 75°.
 - Built-in front safety tread footstep.
 - Sliding legrest extension with removable pad attached with VELCRO® brand fasteners.
 - Heavy duty concealed stirrups adjust in length and are laterally adjustable in (4) positions. Vinyl coated heelrests.
 - Paper Dispenser and Paper Cutter.
 - Hospital grade duplex outlet.
- Note: LD models do not have outlet.**
- Black laminate base with leveling legs.

COLORS:

Laminate – Folkstone Gray

Upholstery– 17 standard vinyl colors

(FOLD OUT FRONT COVER)

A. Series 4423, 4424

SPECIAL FEATURES:

- (2) Drawers in front and (3) Drawers on side as shown.

MODEL	LENGTH	WIDTH	HEIGHT
4423	55" - 71"	26"	33"
4423-LD	Same but less Duplex outlet		
4424	Left Hand Table		
4424-LD	Same but less Duplex outlet		

B. Series 4420, 4421

SPECIAL FEATURES:

- (2) Drawers in front and (1) Door on side with shelf as shown.

MODEL	LENGTH	WIDTH	HEIGHT
4420	55" - 71"	26"	33"
4420-LD	Same but less Duplex outlet		
4421	Left Hand Table		
4421-LD	Same but less Duplex outlet		

OPTIONS

FOR SERIES 4420, 4421, 4423, 4424

- *15 – (1) Safety Rail (Order [2] if both sides)
- *091 – Pelvic Tilt (Controls on both sides)
- *094 – Drawer Warmer
- *097 – Stainless Steel Pan & Cutout

*FACTORY INSTALLED

SEE PAGE 4



Full-comfort pad with VELCRO® brand fasteners detaches easily for cleaning or replacement.



A. Model 4423



B. Model 4420

Shown with option: 097 – Stainless Steel Pan

Metal Paper Dispenser

Upholstered, removable top pad

Pneumatic Air Spring Backrest

Slide-Out leg rest with pad, 21" wide by 16" deep

Storage cabinet with shelf

4-position laterally adjustable stirrups, concealed when not in use

Stainless steel irrigation pan (optional)

Hospital grade duplex outlet (not supplied on "LD" models)

Built-In Safety-Tread Foot Step, 21" wide x 13" deep x 7¾" high