

## A. Model 4710 All-Purpose Hi-Lo Treatment Table

### FEATURES:

- Electric motorized height adjustment from 18½" (wheelchair height) to 35" high.
- "Hands-free" foot control operation.
- Adjustable 26" long headrest has dual gas spring operation and single lever activation.
- Headrest raises to 83°. Lowers to 40° below flat.
- Base has 4½" clearance for most patient hoists.
- Steel frame with powder coated cream finish.
- Stocking upholstery vinyl color is Dove Gray.
- Retractable Locking Casters with lever activation.
- Paper Dispenser at foot end.
- 110V, 60 Hz, 3.15 Amps.



MODEL	LENGTH	WIDTH	HEIGHT
4710	78"	28"	18½" - 35"



A. Model 4710



Model 4710 shown flat



B. Model 4766

*New Design!*

## B. Model 4766 Powermatic® Procedure Table with Gas-Spring Backrest



- 2" High-Density urethane foam vinyl top.
- "Hands Free" foot control adjust height from 26" to 38".
- Includes Paper Dispenser & Cutter.
- Hospital Grade Plug. 110V, 60Hz, 3.5Amps. U.L. Listed.
- New contoured base plate with (4) levelers and satin aluminum trim.
- Gas-spring pneumatic 30" backrest raises to 75° with (1) lever control.

MODEL	LENGTH	WIDTH	HEIGHT
4766	76"	27"	26" - 38"

## Model 4767 (not shown) Powermatic® Procedure Table with Power Backrest



- Same as Model 4766 except Power 30" backrest raises to 75° with foot control.

MODEL	LENGTH	WIDTH	HEIGHT
4767	76"	27"	26" - 38"

## C. Model 4707 Electric Hi-Lo Treatment Table



- Adjustable height motorized hydraulic cylinders are built into each (4) corner steel leg.
- "Hands Free" foot control activation.
- Apron and shelf are almond laminate.
- Black powder coated legs with levelers.
- 110V, 60Hz, 2.4 Amps.

Upholstery – 17 standard vinyl colors.

MODEL	LENGTH	WIDTH	HEIGHT
4707	72"	30"	27"-37"

## Model 4705 (not shown) Crank Hi-Lo Treatment Table



- Same as #4707 except hand crank operation with Fold-away handle.

MODEL	LENGTH	WIDTH	HEIGHT
4705	72"	30"	27"-37"

### OPTIONS FOR MODELS 4705, 4707

- \*06 – Backrest, adjustable
- \*07 – Nose Cutout
- \*27 – Paper Dispenser & Cutter Combo
- \*99L – Left Hand Table Model



C. Model 4707

Shown in left hand view with options:  
06 – Backrest  
27 – Paper Cutter & Dispenser