

Common Features

HI-LO TABLES – MODELS 4701, 4702, 4703, 4704:

- Hydraulic cylinders are built into each (4) Corner Telescoping Legs.
- Black powder coated Legs with Levelers for uneven floors.
- Vinyl upholstered tops with 2" High-Density urethane foam.
- Includes Paper Dispenser and Cutter.
- 16"L legrest extends overall top length to 68".



COLORS: (FOLD OUT FRONT COVER)

Upholstery – 9 Pro-Form™ vinyl colors.



A. Model 4703-722

A. Model 4703 Crank Hydraulic Changing/ Treatment Table



SPECIAL FEATURES:

- Crank hydraulic adjustment from 20" wheelchair height to 30" high.

MODEL	LENGTH	WIDTH	HEIGHT
4703	52" - 68"	30"	20" - 30"
4703-722	52" - 68"	30"	20" - 30"

B. Model 4702 Electric Hi-Lo Changing/ Treatment Table



SPECIAL FEATURES:

- Electrically activated hydraulic cylinders are built into each (4) Corner Telescoping Legs.
- Motorized table lowers to 20" wheelchair height. Elevates up to 30" high.
- "Hands-Free" foot control operation.

MODEL	LENGTH	WIDTH	HEIGHT
4702	52"-68"	30"	20" - 30"



B. Model 4702

C. Model 4701 Crank Hydraulic Changing/ Treatment Table with Drawers



SPECIAL FEATURES:

- Crank hydraulic adjustment from 27"H to 37"H.
- Hydraulic cylinders are built into each (4) corner telescoping legs.

MODEL	LENGTH	WIDTH	HEIGHT
4701	52"-68"	30"	27" - 37"

Model 4704 (not shown) Electric Hi-Lo Changing/ Treatment Table w/ Drawers



SPECIAL FEATURES:

- Electrically activated hydraulic cylinders are built into each (4) Corner Telescoping Legs.
- Motorized table lowers to 27". Elevates up to 37"H.
- "Hands-Free" foot control operation.

MODEL	LENGTH	WIDTH	HEIGHT
4704 (not shown)	52"-68"	30"	27" - 37"



C. Model 4701