

# REHAB SHOWER COMMODES

Traditionally, rehab shower commodes have been constructed of steel or stainless steel. These two materials had their drawbacks. Steel rusts, stainless steel is expensive and both are heavy. Ezee Life™ commodes are constructed of high quality aluminum that won't rust, is lighter than steel and lower in cost than stainless steel.



New locking flip up arms.



Optional 4 Way Seat



Anti tippers  
on all models



Height adjustable  
on all models



Open Front Seat & Insert



Standard on all sizes

Closed Front Seat & Top



Available on 180, 185 & 186 Models

4-Way Open Seat



Available on 180

Removable snap on seat (all models)



## STANDARD FEATURES:

- Lightweight aluminum **NO RUST**
- Footrest, adjustable height & removable
- Removable insert cushion for showering or toileting
- Quick Release axles (self propelled models)
- Four 5" locking casters
- Slide on pail, removable
- Flip back arms
- Maintenance free tires
- Calf strap included
- Molded handrims
- Super soft seat
- Easier rolling
- Seat belt
- All adjustments are tool-less
- 1 Year full warranty
- 2 Year warranty on frame

# TILT REHAB SHOWER COMMODES

Traditionally, rehab shower commodes have been constructed of steel or stainless steel. These two materials had their drawbacks. Steel rusts, stainless steel is expensive and both are heavy. Ezee Life™ commodes are constructed of high quality aluminum that won't rust, is lighter than steel and lower in cost than stainless steel.



New locking flip up arms.



Optional 4 Way Seat



Anti tippers  
on all models



Height adjustable  
on all models



Removable snap on seat (all models)



Open Front Seat & Insert



Standard on all sizes

Closed Front Seat & Top



Available on 190 & 195 Commodes

4-Way Open Seat



Available on 190

## STANDARD FEATURES:

- Lightweight aluminum **NO RUST**
- Footrest, adjustable height & removable
- Removable insert cushion for showering or toileting
- Quick Release axles (self propelled models)
- Adjustable headrest up/down/forward/back
- Four 5" locking casters
- Slide on pail, removable
- Flip back arms
- Maintenance free tires
- Calf strap included
- Molded handrims

- Super soft seat
- Easier rolling
- Seat belt
- All adjustments are tool-less
- 1 Year full warranty
- 2 Year warranty on frame

COMMODOE MODELS	
Model Codes	Seat Width
150	15"
180	18"
185	20"
186	22"
Model numbers indicate the seat size of the commode.	

OPTIONAL EQUIPMENT AND FEATURES		
Option Codes	Option Name	Availability
4W	4 Way seat	All 18" Commodes
CSF18	18" Closed Front Seat	All 18" Commodes
CSF 20	20" Closed Front Seat	All 20" Commodes
CSF 22	22" Closed Front Seat	All 22" Commodes
22	22" Rear Wheels	15" Commodes
24	24" Rear Wheels	18", 20" & 22" Commodes
LSC	Long Stem Casters	All Commodes
*Add the option codes above to the model numbers to select		
*Open Front Seat and 5" Rear Casters Standard on All Models		

SPECIFICATIONS	150, 150-22	180, 180-4W, 180-24, 180-4W-24	185 Wide, 185-24	186 X-Wide, 186-24
Frame Material	Aluminum	Aluminum	Aluminum	Aluminum
Weight Capacity	286 lbs (130 kgs)	286 lbs (130 kgs)	385 lbs (175 kgs)	440 lbs (200 kgs)
Rear Wheels	5" locking casters	5" locking casters	5" locking casters	5" locking casters
Optional Rear Wheels	22" (150-22)	24" (180-24)	24" (185-24)	24" (186-24)
Seat Width	15"	18"	20"	22"
Seat Depth	15"	17"	17"	17"
Foot Rest Hanger Angle	60°	60°	60°	60°
Foot Rest Length	14" to 20"	14" to 20"	14" to 20"	14" to 20"
Back Height	18"	18"	18"	18"
Product Weight	26.4 lbs (12 kgs)	28.4 lbs (12.9 kgs)	29 lbs (13.2 kgs)	32.3 lbs (14.7 kgs)
Arm Height	11" (from seat)	11" (from seat)	11" (from seat)	11" (from seat)
Arm Length	13"	13"	13"	13"
Arm Pad Length	9"	9"	9"	9"
Tilt Range Posterior	n/a	n/a	n/a	n/a
Rear Internal Width	12"	15"	17"	19"
Seat Heights (Standard casters)	19.3", 20.3", 21.3"	19.3", 20.3", 21.3"	19.3", 20.3", 21.3"	19.3", 20.3", 21.3"
Seat Heights (Long stem casters)	21", 22", 23"	21", 22", 23"	21", 22", 23"	21", 22", 23"
Seat Clearance Height <sup>1</sup> (Standard caster stems)	16.3", 17.3" 18.3"	16.3", 17.3" 18.3"	16.3", 17.3" 18.3"	16.3", 17.3" 18.3"
Seat Clearance Height <sup>2</sup> (Long caster stems)	18", 19", 20"	18", 19", 20"	18", 19", 20"	18", 19", 20"
DIMENSIONS WITH 5" REAR CASTERS	150	180	185 Wide	186 X-Wide
Overall Width	19"	22"	24"	26"
Overall Length <sup>1</sup>	37"	37"	37"	37"
DIMENSIONS WITH 24" REAR WHEELS	150	180	185 Wide	186 X-Wide
Overall Width	23"	26"	28"	30"
Overall Length <sup>1</sup>	42"	42"	42"	42"

1) Length with all parts extended (wheels highest height, anti-tip bars fully extended, footrests in place)

2) Seat clearance is the measurement from the bottom of the pail rail to the floor.

COMMODE MODELS	
Model Codes	Seat Width
155	15"
190	18"
195	20"
Model numbers indicate the seat size of the commode.	

OPTIONAL EQUIPMENT AND FEATURES		
Option Codes	Option Name	Availability
4W	4 Way seat	All 18" Commodes
CSF18	18" Closed Front Seat	All 18" Commodes
CSF20	20" Closed Front Seat	All 20" Commodes
22	22" Rear Wheels	15" Commodes
24	24" Rear Wheels	18" & 20" Commodes
LSC	Long Stem Casters	All Commodes
*Add the option codes above to the model numbers to select		
*Open Front Seat and 5" Rear Casters Standard on All Models		

SPECIFICATIONS	155 Tilt, 155-22	190 Tilt, 190-4W, 190-24, 190-4W-24	195 Wide Tilt, 195-24
Frame Material	Aluminum	Aluminum	Aluminum
Weight Capacity	286 lbs (130 kgs)	300 lbs (136 kgs)	325 lbs (148 kgs)
Rear Wheels	5" locking casters	5" locking casters	5" locking casters
Optional Rear Wheels	22" (155-22)	24" (190-24)	24" (195-24)
Seat Width	15"	18"	20"
Seat Depth	15"	17"	17"
Foot Rest Hanger Angle	60°	60°	60°
Foot Rest Length	14" to 20"	14" to 20"	14" to 20"
Back Height	18"	18"	18"
Product Weight	32.5 lbs (14.8 kgs)	36 lbs (16.4 kgs)	38 lbs (17.3 kgs)
Arm Height	11" (from seat)	11" (from seat)	11" (from seat)
Arm Length	13"	13"	13"
Arm Pad Length	9"	9"	9"
Tilt Range Posterior	20°	20°	20°
Rear Internal Width	15"	15"	17"
Seat Heights (Standard casters)	21", 22", 23"	21", 22", 23"	21", 22", 23"
Seat Heights (Long stem casters)	23", 24", 25"	23", 24", 25"	23", 24", 25"
Seat Clearance Height <sup>2</sup> (Standard caster stems)	18", 19", 20"	18", 19", 20"	18", 19", 20"
Seat Clearance Height <sup>2</sup> (Long caster stems)	20", 21", 22"	20", 21", 22"	20", 21", 22"
DIMENSIONS WITH 5" REAR CASTERS	155 Tilt	190 Tilt	195 Wide Tilt
Overall Width	19"	22"	24"
Overall Length <sup>1</sup>	40"	40"	40"
DIMENSIONS WITH 24" REAR WHEELS	155 Tilt	190 Tilt	195 Wide Tilt
Overall Width	23"	26"	28"
Overall Length <sup>1</sup>	45"	45"	45"

1) Length with all parts extended (wheels highest height, anti-tip bars fully extended, footrests in place)

2) Seat clearance is the measurement from the bottom of the pail rail to the floor.

Pictures and specifications may differ slightly.

## SEAT OPTIONS FOR THE EZEE LIFE REHAB SHOWER COMMODES

EZee Life™ now offers up to three different seat options for its aluminum rehab shower commodes, depending on the model, to meet the needs of a large range of users and caregivers.

Our commodes are fitted with a carefully chosen premium quality padded seat, back and arms to enhance client comfort and postural support while being able to withstand the rigours of use and vigorous cleaning solvents found in institutions and homes.

All our seats have a unique tool-less design that can be easily removed from the frame for cleaning while still providing durability and stability. The ease of cleaning these rehab commodes make them one of the most effective commodes available for situations where infection control is paramount.

### OPEN FRONT SEAT

The open front commode seat is the standard seat found on all EZee Life™ aluminum rehab shower commodes. This seat allows hygiene access by the user or a caregiver from the front of the commode and has a padded insert that fills in the commode opening for those who don't want or need to use the commode feature.

The open front seat is available on all models and sizes.



### CLOSED FRONT SEAT


The closed front commode seat is an optional seat that fits several models of EZee Life™ aluminum rehab shower commodes. The closed front is for users who do not require hygiene access while seated. This seat comes with a full padded top cover for users who don't need to have a commode opening all the time but want the option to switch back and forth.

The closed front seat is available on 180, 185, 190 and 195.



### 4-WAY SEAT

The 4-Way Seat is a standard open front commode seat that can be mounted onto the commodes facing the front, rear, left or right side. This seat option allows hygiene access from any direction. The seat can be changed from client to client and from caregiver to caregiver depending on preference and needs.



The 4-Way Seat is only available on the 180 and 190 commode models.



4-WAY FRONT



4-WAY LEFT



4-WAY REAR



4-WAY RIGHT