

US Retail Price List

Account Information

Date	Account Number
Account Name	Contact Name
Phone	Fax
Email Address	
PO Number	Mark For

Quote

Order

Ship To Information

Name		
Address		
City	State	Zip
Phone		

# Hoggi® Bingo Evolution



## Standard Features

- 45 Deg. Tilt in Space
- Front or Rear Facing
- Height Adj. Push Handle
- Cargo Basket
- Reclining Back to 160 Deg.
- Fully Horizontal Back (180 deg) w / combo of Tilt & Recline
- Angle Adj. Legrest
- Compact Folding Design
- Adjustable Seat and Back Shell

## Specifications

### Mobility Base

	Size 1	Size 2	Size 2XL
Push Handle Height	30-44 in.	30-44 in.	27-46.5 in.
Overall Width	25 in.	26.5 in.	26.5 in.
Overall Length	36 in.	36 in.	36 in.
Wheel Diameter, Front	7.5 in.	7.5 in.	7.5 in.
Wheel Diameter, Rear	11.4 in.	11.4 in.	11.4 in.
Tilt in Space (degrees)	45 deg	45 deg	45 deg
Folded Size with Wheels	31.5x25x15 in.	31.5x27x15 in.	31.5x27.15 in.
Folded Size with Wheels Removed	31.5x25x12 in.	31.5x26.5x12 in.	31.5x26.5x12 in.
Weight	19 lbs.	19 lbs.	19 lbs.
Total Load Capacity of Mobility Base (Seat + Rider)	132 lbs.	132 lbs.	132 lbs
** Complex Respiratory Base Total Weight Capacity: 154 lbs.			
Turning Radius	53 in.	55 in.	55 in.

### IUS Seating System

Seat Depth	6.3 to 13.4 in.	8.6 to 15.75 in.	10.6 to 17.7 in.
Seat Width	7 to 11.8 in.	9 to 13.8 in.	9 to 13.8 in.
Lower Leg Length	7.5 to 14.6 in.	7.5 to 15.75 in.	7.5 to 15.75 in.
Backrest Height	17.7 to 25 in.	21 to 29.5 in.	24.4 to 32.7 in.
Height of Slots for Shoulder Straps	11 to 18 in.	13 to 21.5 in.	16.5 to 24.8 in.
Recline	90° to 160°	90° to 160°	90° to 130°
**Back completely horizontal (180 deg) with combo tilt & recline			
Weight	14.1 lbs.	16.3 lbs	18.3 lbs
Load Capacity of Seat Unit (Rider Weight Limit)	110 lbs.	110 lbs.	110 lbs.

Transit Brackets Crash tested according to ISO 7176-19 ANSI/RESNA WC/Vol. 1. – Section 19



Size 1	Mobility Base	E1234	3232-1000	<input type="checkbox"/>	\$2,540	All Items Required
	Reclining Backrest	E1014	3231-0155	<input type="checkbox"/>	\$520	
	Elevating Legrest	E0990	3231-2400	<input type="checkbox"/>	\$285	
	Bingo Evolution IUS Seat Cushion	E2292	3231-5100-1	<input type="checkbox"/>	\$685	
	Bingo Evolution IUS Ht. Adjustable Back Cushion	E2293	3231-5100-2	<input type="checkbox"/>	\$685	

Size 2	Mobility Base	E1234	3234-1000	<input type="checkbox"/>	\$2,540	All Items Required
	Reclining Backrest	E1014	3233-0155	<input type="checkbox"/>	\$520	
	Elevating Legrest	E0990	3233-2400	<input type="checkbox"/>	\$285	
	Bingo Evolution IUS Seat Cushion	E2292	3233-5100-1	<input type="checkbox"/>	\$685	
	Bingo Evolution IUS Ht. Adjustable Back Cushion	E2293	3233-5100-2	<input type="checkbox"/>	\$685	

Size 2XL	Mobility Base	E1234	3235-1000	<input type="checkbox"/>	\$2,860	All Items Required
	Reclining Backrest	E1014	3235-0155	<input type="checkbox"/>	\$520	
	Elevating Legrest	E0990	3233-2400	<input type="checkbox"/>	\$285	
	Bingo Evolution IUS Seat Cushion	E2292	3235-5100-1	<input type="checkbox"/>	\$795	
	Bingo Evolution IUS Ht. Adjustable Back Cushion	E2293	3235-5100-2	<input type="checkbox"/>	\$795	
	2XL Extended Push Handle	K0108	SA-0042-MP	<input type="checkbox"/>	\$228	

Transit Package (Required for all bases)	K0108	3231-7300	<input checked="" type="checkbox"/>	\$335
------------------------------------------	-------	-----------	-------------------------------------	-------

Frame Color: Iron Sparkle 024  Matte Black 032

Upholstery Color: Red 345  City Stone 338  Tahoe Blue 339  Hot Pink 340  \$4

**One** selection below is **Required**.

Select **either** (A) Hip Guides or (B) Armrests:

		Size 1			Size 2 & 2XL			
<b>A</b>	Hip Guides (pair) (requires hardware)	E0956	3231-2600-323	<input type="checkbox"/>	\$105	Size 2 3233-2600-323	<input type="checkbox"/>	\$115
						Size 2XL 3235-2600-323	<input type="checkbox"/>	\$125
	Adjustable Removable Hardware (pair)	E1028	HHWDR	<input type="checkbox"/>	\$300	HHWDR	<input type="checkbox"/>	\$300
<b>B</b>	Height Adjustable Armrest	E0973	3231-3201	<input type="checkbox"/>	\$434	3233-3201	<input type="checkbox"/>	\$434
	**Not compatible with Hip Supports listed above Hip Pads for inside armrest plates (Gray, optional)	E0956	3201-5050-323	<input type="checkbox"/>	\$196	3203-5050-323	<input type="checkbox"/>	\$196




				Size 1			Size 2 & 2XL			
Upper Body Accessories	Headrest Lateral Support Pads (gray)		E0955	3231-7600-323	<input type="checkbox"/>	\$328	3233-7600-323	<input type="checkbox"/>	\$328	
	Lateral Trunk Supports (pair)		E0956	Choose	Sz 1 3231-7700-323	<input type="checkbox"/>	\$165	Sz 2 3233-7700-323	<input type="checkbox"/>	\$175
					Sz 2 3233-7700-323	<input type="checkbox"/>	\$175		Sz 3 3234-7700-323	<input type="checkbox"/>
	Adjustable Hardware (pair) **required		E1028		LHWDR	<input type="checkbox"/>	\$330	LHWDR	<input type="checkbox"/>	\$330
	Chest Harness (Gray)		E0960	Choose	Sm 3201-5600-323	<input type="checkbox"/>	\$263	Med 3201-5700-323	<input type="checkbox"/>	\$263
			Med 3201-5700-323		<input type="checkbox"/>	\$263	Lg 3203-5600-323		<input type="checkbox"/>	\$274
5 Point Harness		E0960		3201-6200	<input type="checkbox"/>	\$259	3201-6200	<input type="checkbox"/>	\$259	

Lower Body	Medial Thigh Support (Gray)		E0957	Choose	Sm 3231-7800-323	<input type="checkbox"/>	\$140	Med 3231-7900-323	<input type="checkbox"/>	\$140	
					Med 3231-7900-323	<input type="checkbox"/>	\$140		Lg 3233-7800-323	<input type="checkbox"/>	\$140
	Swing Away Removable Hardware **required		E1028		MTHWDR	<input type="checkbox"/>	\$280	MTHWDR	<input type="checkbox"/>	\$280	
	Groin Strap/Leg Harness		E0978		3231-5300	<input type="checkbox"/>	\$239	3233-5300	<input type="checkbox"/>	\$249	
	Lap Belt		E0978		3201-6100	<input type="checkbox"/>	\$171	3201-6100	<input type="checkbox"/>	\$171	
	4 Point Lap Belt		E0978		3231-6702	<input type="checkbox"/>	\$269	3231-6703	<input type="checkbox"/>	\$269	
Footstraps (pair)		E0952		3201-6300	<input type="checkbox"/>	\$72	3201-6300	<input type="checkbox"/>	\$72		
Ankle Cuffs (pair)		K0108		3201-6400	<input type="checkbox"/>	\$267	3203-6400	<input type="checkbox"/>	\$277		






				Size 1			Size 2 & 2XL			
Additional Options	Canopy with Rain cover		K0108	3201-9300	<input type="checkbox"/>	\$425	3201-9300	<input type="checkbox"/>	\$425	
	Tray		E0950	Gray	<input type="checkbox"/>	\$251	Gray	<input type="checkbox"/>	\$251	
	**Requires hardware below									
	Tray Mounting Hardware, Removable		E1028	Choose	3231-8201	<input type="checkbox"/>	\$262	3233-8201	<input type="checkbox"/>	\$262
	Tray Mounting Hardware, Removable & Height Adjustable		E1028		3231-8203	<input type="checkbox"/>	\$336	3233-8203	<input type="checkbox"/>	\$336
	Grab rail with upholstery		K0108		3231-8100	<input type="checkbox"/>	\$171	3233-8100	<input type="checkbox"/>	\$182
	4" Footrest Risers		K0108		riser4	<input type="checkbox"/>	\$78	riser4	<input type="checkbox"/>	\$78
	Double female threaded, solid aluminum posts									
	6" Footrest Risers		K0108		riser6	<input type="checkbox"/>	\$89	riser6	<input type="checkbox"/>	\$89
Extended Length Elevating Legrest Shaft (3" Longer)		K0108		SA-0300/2-MP	<input type="checkbox"/>	\$191	SA-0300/2-MP	<input type="checkbox"/>	\$191	
Calf Plate Upholstery Insert (Metal plate and foam inside upholstery)		K0108		3232-5190	<input type="checkbox"/>	\$107	3234-5190	<input type="checkbox"/>	\$107	

Cobra Indoor Base		K0108	3004-0000-024	<input type="checkbox"/>	\$2,747	3005-0000-024	<input type="checkbox"/>	\$2,747
			Iron Sparkle			Iron Sparkle		
			3004-0000-032	<input type="checkbox"/>	\$2,747	3005-0000-032	<input type="checkbox"/>	\$2,747
			Matte Black			Matte Black		

Continued on next page



			Size 1			Size 2 & 2XL		
Respiratory & Complex Care Equipment Options	Oxygen Tank Holder † 	E2208	3231-4202	<input type="checkbox"/>	\$636	3231-4202	<input type="checkbox"/>	\$636
	IV Pole † (Removable and height adjustable)	K0105	7000-0081	<input type="checkbox"/>	\$275	7000-0081	<input type="checkbox"/>	\$275
	4" Seat Riser Modification 	K0108	SA-0298-KP	<input type="checkbox"/>	\$687	SA-0298-KP	<input type="checkbox"/>	\$687
	**This modification raises the entire seat by 4". This allows greater storage capacity for medical equipment and eases transfers and caregiving. Can be selected with or without complex respiratory base modification.							
	Respiratory Equipment Tray † 	E1029	US3239-vent	<input type="checkbox"/>	\$495	US3239-vent	<input type="checkbox"/>	\$495
**Must Order 3239- 10 06- 024. The vent tray consists of Three separate panels: Front = 13"w x 8"d Middle = 19.5"w x 11.25"d, Rear = 19.75"w x 11.3"d Includes mounting tower for LTV***								
Complex Respiratory Base Modification **This modification includes a longer push handle, modified frame. This is a modified Size 2 Mobility Base that will fit both Size 1 and Size 2 seat shells. Seat adapter required when Size 1 seat selected. Increased weight capacity of chassis to 154lbs!	K0108	3239-1006	<input type="checkbox"/>	\$1650	3239-1006	<input type="checkbox"/>	\$1650	
Wheel Braking System & Lock	E2228	SA-0219-KP Adapter Required with Size 1 Seat	<input type="checkbox"/>	\$432				
Right and Left <u>required</u> w/ventilator tray. Std brakes not available.		WBSL Right	<input type="checkbox"/>	\$1100	WBSL Right	<input type="checkbox"/>	\$1100	
		WBSL Left	<input type="checkbox"/>	\$1100	WBSL Left	<input type="checkbox"/>	\$1100	

† Oxygen Tank Holder, IV Pole and Respiratory Equipment Tray have not been crash tested  
 \*\*\*Ventilator specific hardware and adapter plates for mounting tower not included

