

HOGGI
MADE IN GERMANY

BINGO EVOLUTION

REHAB PUSHCHAIR

Sizes 1-2 • 2 XL • 3
Up to 110 lbs


RUGGED &
ADAPTABLE
FOR ACTIVE
FAMILIES



	Size 1	Size 2	Size 2XL	Size 3
MOBILITY BASE	3232-1000	3234-1000	3235-1000	
Push Handle Height	30" - 43"	30" - 45"	34" - 49"	34" - 49"
Overall Width/Length	25"/36"	27"/36"	27"/36"	29"/38"
Wheel Diameter Front/Rear	7.5"/11"	7.5"/11"	7.5"/11"	7.5"/11"
Tilt-in-Space	45°	45°	45°	45°
Folded size w/Wheels	32" x 25" x 15"	32" x 27" x 15"	32" x 27" x 15"	32" x 29" x 17"
Base only Weight	19 lbs	19 lbs	19 lbs	20 lbs
Rider Weight Limit	110 lbs	110 lbs	110 lbs	110 lbs
Total Weight Capacity (Seat + Rider) Standard Base	132 lbs	132 lbs	132 lbs	132 lbs
Total Weight Capacity (Seat + Rider) Respiratory Base	154 lbs	154 lbs	154 lbs	154 lbs
IUS SEATING SYSTEM				
Seat Depth	6" - 13"	9" - 16"	11" - 18"	11" - 18"
Seat Width	7" - 12"	9" - 14"	9" - 14"	14" - 17"
Lower Leg Length* (without optional risers)	8" - 15"	8" - 16"	8" - 16"	8" - 18"
Backrest Height	18" - 25"	21" - 30"	24" - 33"	24" - 33"
Height of Slots for Shoulder Straps	11" - 18"	13" - 22"	17" - 25"	17" - 25"
Recline 90° - 160° *with recline, backrest reaches 180°		90° - 130°		

	Size 1	Size 2	Size 3
OT PLUS SEAT SHELL			
*Separate Back, Seat and accessories supplied by dealer.			
Seat Depth	6" - 13"	9" - 16"	11" - 18"
Seat Width	Custom - 12"	Custom - 14"	Custom - 17"
* When ordering afermarket seating, max solid base width	11.75"	13.75"	16.75"
Back Cane Tubing Diameter	1"	1"	1"
Back Cane Height (Pan to top of cane)	15.5"	17"	23"
Lower Leg Length* (without optional risers)	8" - 15"	8" - 16"	8" - 18"
Recline - 90° - 160° *with recline, backrest reaches 180°			90° - 130°
Seat only Weight	14 lbs	16 lbs	17 lbs

Transit Brackets Crash tested according to IS) 7176-19 ANSU/RESNA WC / VOL. 1 - Section 19



distributors of
HOGGI



distributors of
HOGGI



Berry Soft



Black



Turquoise Soft



Gray Soft

SUPPORTS

Headrest
Truck Support
Hip & Thigh

SEAT ANGLES

45° Tilt
160° Reclining back

SHOCKS

Durable
Tool-Free
4-Settings

ADJUSTABLE LEG REST

Flip-up footplate

CARGO BASKET

Standard Feature

PUSH BAR

Adjustable Height

TRANSPORTATION BRACKETS

Required Component

LARGE WHEELS

Flat-Free
Stylish

COBRA BASE



Indoor Hi-Low


Full Tilt/Recline



Headrest Pads



Highrise Headrest Pads



Anatomic Head Support



Trunk Supports



Chest-shoulder Harness




5-point Harness



Lap Belt



4-point Lap Belt



Grab Rail With Upholstery



Hip Supports
(required component)



Removable Abduction Block



Groin Strap



Armrests



Therapy Tray



Ankle Cuffs



Foot Straps



Canopy with Rain Cover



3 Piece Removable Vent Tray



Oxygen Bottle Holder



Frame Riser

Bingo Evolution OT Plus Frame for After Market Seating

The OT Plus seat shell has all the same Tilt and Recline capabilities but allows the attachment of after market seating of your choice to the 1" round tubing.

(Seating depicted to the right is for example purposes only and purchased by dealer from a separate manufacturer)

