



### FEATURES

- Accommodate users weighing 20 - 102 lb and up to 60" tall
- Durable, machine washable, comfortable padded cover in three color options
- Crash tested to FMVSS #213
- Certified for use in aircraft
- Positioning options available
- Optional seat depth extender allows car seat to grow with the child

### INCLUDES

- Soft-padded headwings
- Memory foam back/seat padding reduces sores and fatigue
- Standard 5-point harness restraint for better protection than vehicle lap/shoulder belt
- Crotch strap pad
- Lateral positioning pads
- Dual-hook tether strap

### UPHOLSTERY



Convertible Pink  
(RR)



Roadster Red  
(RR)



Speedway Gray  
(SC)

### ACCESSORIES



Seat Depth Extender



Buckle Guard



Retainer Clip Guard



Extensor Thrust  
Wedge

## ORDER FORM

Dealer: \_\_\_\_\_ Marked For: \_\_\_\_\_ Account #: \_\_\_\_\_

Date: \_\_\_\_\_ PO#: \_\_\_\_\_ Ordered By: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_ E-mail: \_\_\_\_\_

Ship to address: \_\_\_\_\_








Must select a package below. Cannot configure your own.

BASE UNIT	CSI-2000W	CSI-2000WA	CSI-2000WM	CSI-2000WMA	CSI-2000WL	CSI-2000WLA
	\$1,315	\$1,749	\$1,360	\$1,961	\$1,375	\$1,770
Weight Capacity	20 -102 lb	20 -102 lb	20 -102 lb	20 -102 lb	20 -102 lb	20 -102 lb
Height Capacity	60"	60"	60"	60"	60"	60"
Shoulder Harness Heights	11.5", 14.5", 17.5", 20.5"	11.5", 14.5", 17.5", 20.5"	11.5", 14.5", 17.5", 20.5"	11.5", 14.5", 17.5", 20.5"	11.5", 14.5", 17.5", 20.5"	11.5", 14.5", 17.5", 20.5"
Back Height	29"	29"	29"	29"	29"	29"
Seat Depth	10"	10"	12.5"	12.5"	15"	15"
Inside Seat Width at Front	15.5"	15.5"	15.5"	15.5"	15.5"	15.5"
Inside Seat Width at Rear	13"	13"	13"	13"	13"	13"
Overall Width	19"	19"	19"	19"	19"	19"
Overall Depth	11"	11"	11"	11"	11"	11"
Seat-to-Back Angle	88°	88°	88°	88°	88°	88°
Product Weight	11 lb	11 lb	11 lb	11 lb	11 lb	11 lb

### PACKAGE INCLUDES

Seat Depth Extender			CSI-K2025	CSI-K2025	CSI-K2026	CSI-K2026
Abductor		CSI-K2015		CSI-K2018		CSI-K2019

### UPHOLSTERY

Convertible Pink with Gray Mesh						
Roadster Red with Black Mesh						
Speedway Gray with Gray Mesh						

### ACCESSORIES

Lateral Pads (set of 4)	CSI-A2010	\$66
Bolt-in Hardware Kit with Strap (for vehicle manufactured before 2002)	CSI-A2031	\$99
Block Style Trunk Support Pads	CSI-A2008-BK	\$40
Buckle Guard	ZZ-A2075	\$45
Extensor Thrust Wedge, 2"	CSI-A2050	\$138
Retainer Clip Guard	CSZ-A2085	\$50