


Contour™ Wrap-Around Back Support

BATHING



ACCESSORIES



Pelvic Support Belt



The Contour™ Wrap-Around Bath Support provides mild positioning assistance and allows the user to sit on the floor of the bathtub while bathing.

BENEFITS

Adapts to User Needs

FEATURES

- Multiple sizes, Lo-Back™ and Hi-Back™ padded support options available with harness
- Orthopedic back height easily adjusts vertically

Durable & Long-Lasting

- Rustproof frame and tear-resistant mesh fabric
- Bacterial and fungal-resistant ComfortMesh™ seat
- Latex and silicone-free

Built-In Growth

- Mounting Base allows upgrade or replacement of support back

Versatile Design

- Wide-stance stable base fits in any tub
- Compact for easy storage or transport

For more information, call us today!

Contour™ Wrap-Around Bath Support



Dealer: _____ Marked For: _____ Account #: _____
 Date: _____ PO #: _____ Ordered By: _____
 Phone: _____ Fax: _____ Email: _____
 Ship to Address: _____

BASE UNIT

LO-BACK WRAP-AROUND

HI-BACK WRAP-AROUND

	BCWL-7210S	BCWL-7310M	BCWH-7610S	BCWH-7710M
	\$683.10	\$745.72	\$796.95	\$802.64
Product Size	Small	Medium	Small	Medium
Weight Capacity	N/A	N/A	N/A	N/A
Overall Depth	29"	29"	29"	29"
Seat Depth	14.5"	14.5"	14.5"	14.5"
Overall Seat Width	18"	18"	18"	18"
Inside Seat Width	15"	15"	15"	15"
Inside Back Width at Top	7"	11"	7"	8"
Inside Back Width at Bottom	7"	11"	12"	12"
Back Height Range	9–14"	10–14"	17–22"	19–24"
Product Weight	5.5 lb.	6.5 lb.	9.5 lb.	9.75 lb.

ACCESSORIES (fit all sizes of bath chairs)

Pelvic Positioning Belt \$79.70 BCWZ-A7830