

# Marlin™ Bath Lift

BATHING



## ACCESSORIES



Head and Trunk  
Support Pads



Contoured Chest  
Harness



Swivel and  
Transfer Aid



Rapid Dry  
Positioning Belt



The compact and easy-to-use Marlin™ Bath Lift features a reclining back and rechargeable battery-operated lift mechanism that is easy to transport or store when not in use.

## BENEFITS

## FEATURES

### User Accommodation

- Seat height adjusts from 2.5–19" and back reclines from 10° off vertical to 50° at its lowest position
- Multiple positioning accessories available for postural support

### Compact & Versatile

- Weighs 20.5 lb. and fits most bathtubs
- Supports up to 300 lb. and includes side flaps for easy lateral transfers

### Easy Operation

- Waterproof, simple-to-use hand control safely and gradually lowers user into bathtub
- Safety mechanism prevents lift from lowering if battery charge is not sufficient to raise it back up

### Travel-Friendly

- Tool-free assembly and installation
- Disassembles for easy storage or transport

For more information, call us today!

©2021 Medical Depot, Inc. dba Drive DeVilbiss Healthcare  
Trademarks owned by Medical Depot

Marlin\_SS\_rev\_04\_22\_2021

Dealer: \_\_\_\_\_ Marked For: \_\_\_\_\_ Account #: \_\_\_\_\_

Date: \_\_\_\_\_ PO #: \_\_\_\_\_ Ordered By: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_ Email: \_\_\_\_\_

Ship to Address: \_\_\_\_\_

## BASE UNIT

\$2,271.31

**BLM-8200**

<b>Weight Capacity</b>	300 lb.
<b>Seat Width</b>	13.8"
<b>Seat Depth</b>	16.5"
<b>Seat Height Range</b>	2.5–19"
<b>Back Height</b>	25.2"
<b>Folded Dimensions</b>	15" W x 28" L x 7.75" H
<b>Battery Specification</b>	12v
<b>Duration of Lifting / Lowering</b>	50 seconds per load
<b>Lift per Full Charge</b>	Approximately 6 lifts with an average 165 lb.
<b>Automatic Shutoff</b>	Maximum 30 minutes
<b>Reclines</b>	10–50°
<b>Product Weight</b>	20.5 lb.

## ACCESSORIES

<b>Head Support Pads</b>	\$176.47	BLM-A8875
<b>Lateral Support Pads</b>	\$176.47	BLM-A8830
<b>Contoured Chest Harness</b>	\$204.93	BLM-A8207
<b>Rapid Dry Positioning Belt</b>	\$239.09	BCZ-A8832M 4" Wide, Fit Range: 22–68"
	\$239.09	BCZ-A8832L 5" Wide, Fit Range: 25–72"
<b>Swivel &amp; Transfer Aid</b>	\$222.01	BLZ-A8202