



FEATURES

- All-in-one welded steel frame available in red and blue
- Plastic seat and contoured back
- Durable snap-on toilet seat
- Height adjustable seat allows child to firmly place their feet on the floor for comfort, safety and security
- Adjustable hook & loop chest strap for additional support
- Molded plastic armrests provide additional comfort
- Rubber non-skid tips ensure that the legs will not mark the floor
- Complete with 7.5 qt removable commode bucket with metal handle and cover



All-in-one welded Steel Frame



Adjustable Hook & Loop Chest Strap



Snap-on Plastic Toilet Seat



Molded Plastic Armrests



Rubber Non-skid Tips



7.5 qt Removable Commode Bucket and Cover

ACCESSORIES



Biodegradable Sanitary Commode Liner

COLORS



Castle Red



Knight Blue

ORDER FORM

Dealer: _____ Marked For: _____ Account #: _____
 Date: _____ PO#: _____ Ordered By: _____
 Phone: _____ Fax: _____ E-mail: _____
 Ship to address: _____

PRODUCT SPECIFICATIONS

	PC 1000 R 1/ea	PC 1000 BL 1/ea
	\$349.00	\$349.00
Seat Depth	15"	15"
Seat Width	13"	13"
Seat Height	10.5" – 14.5"	10.5" – 14.5"
Armrest Height (from floor)	15.5" - 19.5"	15.5" - 19.5"
Armrest Pad	2" (W) x 9" (D)	2" (W) x 9" (D)
Overall Width	17"	17"
Toilet Opening	6.5" (W) x 8" (D)	6.5" (W) x 8" (D)
Weight	9.5 lbs.	9.5 lbs.
Carton Shipping Weight	13 lbs.	13 lbs.
Weight Limit	250 lbs.	250 lbs.
Warranty	Limited Lifetime	Limited Lifetime

ACCESSORIES

Poziflex Chest Strap, Extra Small	PCB-900XS 3" Wide, 12.5" Padded Length	\$85.00
Poziflex Chest Strap, Small	PCB-900S 3.75" Wide, 15.75" Padded Length	\$85.00
Biodegradable Commode Liner	PCB-900S 12/bx, 6 box /cs	\$232.80