

Spirit™ Spica

CAR SEAT



The Inspired by Drive Spirit™ Spica Car Seat features a low and wide profile with positioning wedges to safely transport users with spica casts. The Spirit Spica seat with wipeable cover accommodates a wide range of user heights and weights to meet individual needs.




ACCESSORIES



Positioning Wedge
(2 included)



Memory foam
support



Wipeable,
low maintenance
upholstery



BENEFITS

FEATURES

Built-In Safety & Security

- 5-point adjustable seat harness with included tether strap provides added safety
- Easily installed in most vehicles
- Two seat depth options with a wide weight and height range makes this seat a great choice for multiple clients

Adjusts to User

- Seat depth options of 12" and 14"
- Accommodates users up to 66" tall and has a 25-130 lb weight capacity
- Adjustability and seat configurations can be obtained for different sized clients and Spica casts

Comfort & Convenience

- High quality memory foam on the back and seating surface provides a comfortable environment
- Padded, wipeable Stealth Black upholstery
- 2 wipeable foam seat wedges included standard
- 1 year limited warranty

For more information, call us today!

Spirit™ Spica Car Seat



Dealer: _____ Marked For: _____ Account #: _____
 Date: _____ PO #: _____ Ordered By: _____
 Phone: _____ Fax: _____ Email: _____
 Ship to Address: _____

BASE UNIT

CSS-2400-SPICA
 \$1,599.00

CSS-2401-SPICA
 \$1,650.00

Weight Capacity	25–130 lb	25–130 lb
Max. Height Capacity	66"	66"
Harness Height	11", 13.5", 16", 18.5", 20.5", 23.5" (min. 10")	11", 13.5", 16", 18.5", 20.5", 23.5" (min. 10")
Shoulder Width	18"	18"
Back Height	28"	28"
Seat Depth	12"	16"
Inside Seat Depth Width at Front	17"	17"
Inside Seat Depth Width at Rear	15"	15"
Overall Width	21.5"	21.5"
Seat-to-Back Angle	98°	98°
Product Weight	29 lb.	29 lb.

INCLUDED ACCESSORY

Positioning Wedges: 13.25" x 7.75" x 4.75" (2 included)

