



PL414RBDDA

### FEATURES

- State-of-the art hydraulic reclining mechanism
- Backrest reclines from 90° to 180°
- Double-coated chrome carbon steel frame for maximum durability
- Durable flame retardant nylon upholstery
- Padded flip-back detachable desk armrests
- Composite lightweight 24" Mag-style rear wheels with urethane tires and 6" front rubber tires
- Offset rear wheel position provides additional stability and evenly distributes user's weight in all reclining positions
- Precision sealed front and rear wheel bearings
- Front caster forks are adjustable in three positions
- Standard dual axle provides easy transition of seat to hemi-level
- Push-to-lock wheel locks
- Spreader bar functions as a back stabilizer and a push handle for caregiver
- Standard with height adjustable swing-away elevating leg rests
- Durable and lightweight composite footplates
- Flip-up anti-tippers with wheels



Lightweight 24" Mag-style Rear Wheels




State-of-the art Hydraulic Reclining Mechanism



Backrest Reclines From 90° to 180°

### ACCESSORIES



Contoured Chest Harness



Backpack Style Chest Harness



Pelvic Belt



Ankle Positioners



Calf Panel

## ORDER FORM

Dealer: \_\_\_\_\_ Marked For: \_\_\_\_\_ Account #: \_\_\_\_\_

Date: \_\_\_\_\_ PO#: \_\_\_\_\_ Ordered By: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_ E-mail: \_\_\_\_\_

Ship to address: \_\_\_\_\_

## PRODUCT SPECIFICATIONS

	<b>PL412RBDDA</b>	<b>PL414RBDDA</b>
	\$1,648.92	\$1,648.92
HCPCS	E1236, E1226	E1236, E1226
Riggings	Elevating Leg rest/Foot rest	Elevating Leg rest/foot rest
Arms	Desk Arm	Desk Arm
Seat-To-Floor Height	15"	16"
Back Height	26"	26"
Overall Length w/ Riggings	38"	38"
Overall Height	40.5"	40.5"
Weight Capacity	250 lbs.	250 lbs.
Warranty	Limited Lifetime	Limited Lifetime

## ACCESSORIES

\*Requires Attachment Kit

Contoured Chest Harness, Standard, Front Pull*	PCH-120	XS	S	M	L	XL	\$145.00
Contoured Chest Harness, Standard, Rear Pull*	PCH-220	XS	S	M	L	XL	\$145.00
Backpack Style Harness	PCH-620	XS	S	M	L	XL	\$120.00
<b>Poziform Harness Attachment Kit:</b>							
For Frame Width 11"-15"	PHA-100S						\$155.00
For Frame Width 15"-19"	PHA-100L						\$155.00
2 pt Pelvic Belt with Push Button	PPB-100PB	XS	S	M	L	XL	\$85.00
4 pt Pelvic Belt with Push Button	PPB-500PB	XS	S	M	L	XL	\$110.00
Calf Panel	PCP-280		S	M	L		\$85.00
Ankle Positioner	PAP-100	XS	S	M	L	XL	\$130.00